

Enteral Feeding Tube

Wide Range of Choices to Meet Patient Needs



Please consult package insert for more detailed safety information and instructions for use.

Bard, Ponsky, and Ponsky-Gauderer are trademarks and/or registered trademarks of C. R. Bard, Inc.

All other trademarks are the property of their respective owners.

Copyright © 2016, C. R. Bard, Inc. All Rights Reserved.

BPV/FEED/0616/0012(1)

BARD

Bard Peripheral Vascular, Inc.
1625 W. 3rd Street Tempe, AZ 85281 USA
Tel: 1 480 894 9515 / 1 800 321 4254
Fax: 1 480 966 7062 / 1 800 440 5376
www.bardpv.com

BARD

Percutaneous Endoscopic Gastrostomy (PEG)

BARD® PEG feeding devices offer physicians a wide range of choices to meet specific patient needs

- Soft and compressible internal dome designed to minimize trauma at insertion
- Identification of tube obstructions enhanced with translucent tubing
- Designed for enhanced patient comfort with soft silicone external bolster and tubing



WIDE RANGE OF
choices




PONSKY® PEG Feeding Devices

Radiopaque Dome allows X-ray confirmation

	Product Codes						
	Tube		Dual Port Adaptor	Safety ¹		Non-Safety	
	Outer Diameter	Inner Diameter	Inner Diameter	Pull	Push ^{††}	Pull	Push ^{††}
Deluxe Kits (Snare & Scissors included)	20F	13.8F	10.9F	003029	003030	000792	000793
	28F	19F	13F	003027	-	-	-
Standard Kits (Snare & Scissors NOT included)	20F	13.8F	10.9F	-	-	000330	000331

The design of these enteral feeding products do not restrict their use to enteral feeding devices only. Please check local laws for any applicable usage restrictions.



PONSKY-GAUDERER® PEG Feeding Devices

Radiopaque Dome allows X-ray confirmation
Requires Endoscopic Removal

	Product Codes						
	Tube		Dual Port Adaptor	Safety ¹		Non-Safety	
	Outer Diameter	Inner Diameter	Inner Diameter	Pull	Push ^{††}	Pull	Push ^{††}
Standard Kits (Snare & Scissors NOT included)	20F	13.8F	10.9F	-	-	000329	000328

The design of these enteral feeding products do not restrict their use to enteral feeding devices only. Please check local laws for any applicable usage restrictions.

¹Safety kits include safety engineered sharp devices ^{††}Push=Guidewire

Replacement Devices

Low-Profile Replacement Devices

Non-Balloon

- BARD® Button Gastrostomy Tube and Accessories
- Cosmetically appealing low profile design
- Designed with simple, easy to connect tubing sets
- Obturator facilitates button placement
- Anti-reflux valve helps to provide safer feeding
- Kits include a continuous and bolus feeding tube



F	Stoma Length (cm)	Product Codes			
		BARD® Button Gastrostomy Tube Kit	Decompression Tube	BARD® Button Continuous Feeding Tube (90° adaptor)	BARD® Button Bolus Feeding Tube (straight adaptor)
18F	1.2	000292	000361		
	1.7	000282	000350		
	2.4	000283	000351	000256	000257
	3.4	000284	000352		
24F	1.2	000293	000362		
	1.7	000285	000353		
	2.4	000286	000354	000258	000259
	3.4	000287	000355		
	4.4	000296	000356		
28F	1.5	000261	000357		
	2.7	000262	000358	000268	000269
	4.3	000263	000359		

The design of these enteral feeding products do not restrict their use to enteral feeding devices only. Please check local laws for any applicable usage restrictions.

Replacement Devices



F	Product Codes	
	5cc Balloon	20cc Balloon
12F	000712	-
14F	000714	-
16F	-	000716
18F	-	000718
20F	-	000720
22F	-	000722
24F	-	000724

The design of these enteral feeding products do not restrict their use to enteral feeding devices only. Please check local laws for any applicable usage restrictions.

Standard Profile Replacement Devices Balloon

BARD® Tri-Funnel Replacement Gastrostomy Tube

- Inflatable internal retention balloon
- Radiopaque stripe for X-ray placement verification
- Easy to identify obstructions with translucent silicone tubing
- Meets patient sizing needs with adjustable external bolster
- Separate medication port

Standard Profile Replacement Devices Non-Balloon

PONSKY® Non-Balloon Replacement Gastrostomy Tube

- Allows X-ray placement confirmation with radiopaque dome
- Bold centimeter markings for easy identification
- Easy to identify obstructions with translucent silicone tubing
- Obturator facilitates tube placement
- Meets patient sizing needs with adjustable external bolster
- Separate medication port



Product Code			
Tube		Dual Port Adaptor	PONSKEY® Non-Balloon Replacement Tube
Outer Diameter	Inner Diameter	Inner Diameter	
20F	13.8F	10.9F	000700

The design of these enteral feeding products do not restrict their use to enteral feeding devices only. Please check local laws for any applicable usage restrictions.

Accessories

PONSKEY® Dual Port Feeding Adaptor

- Available for 20F and 28F PONSKEY® devices
- Made of silicone




Product Codes		
PONSKEY® PEG	28F	000318
PONSKEY-GAUDERER® PEG	20F	000333
PONSKEY® PEG	20F	000333

The design of these enteral feeding products do not restrict their use to enteral feeding devices only. Please check local laws for any applicable usage restrictions.

Jejunal Feeding/Gastric Decompression Tube

- Designed for easy insertion with Hydromer™ coating
- Simultaneous feeding and decompression facilitated by two port design



Product Codes			
J-Tube Outer Diameter	Compatible Bard PEG Size	Pull 27 Inch	Push 35 Inch
12F	28F	000733	000734
9F	20F	000319	000732

The design of these enteral feeding products do not restrict their use to enteral feeding devices only. Please check local laws for any applicable usage restrictions.

BARD® PEG Cleaning Brush

- Excellent tube cleaning with nylon brush
- Optimized flushing with cleaning solutions provided through utilization of Luer hub
- Silicone tip helps to protect stomach



Product Code
000396

Stoma Measuring Device

- Expandable/collapsible retention dome facilitates stoma length measurement
- Ease of insertion with consistent tubing diameter



Product Codes	
Used with Button 18F/24F	000842
Used with Button 28F	000843

