

Crutches and replacement parts

After a few days on crutches, you know how uncomfortable they can get. Our answer for sore armpits?

Replacing the traditional foam rubber with **Thermoplastic Rubber (TPR)**—a synthetic blend that lasts long, resists peeling and most important, is soft. Your patients will feel the difference—and so will you, as you help them heal.

Crutches

Cat. no.	Weight capacity	Recommended user height	Description	Packaging
CA901AD	300 lb.	5 ft. 2 in.–5 ft. 10 in. (Adult)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA801ADB	650 lb.	5 ft. 2 in.–5 ft. 10 in. (Adult)	Axillary, heavy duty, push button, adjustable, steel	4 pr/cs
CA901TL	300 lb.	5 ft. 10 in.–6 ft. 6 in. (Tall)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA801TLB	650 lb.	5 ft. 10 in.–6 ft. 6 in. (Tall)	Axillary, heavy duty, push button, adjustable, steel	4 pr/cs
CA901YTH	300 lb.	4 ft. 6 in.–5 ft. 2 in. (Youth)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA901CH	300 lb.	4 ft. 0 in.–4 ft. 6 in. (Child)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA801FYTH	250 lb.	5 ft. 3 in.–6 ft. 3 in. (Tall)	Forearm, push button, adjustable, aluminum	10 pr/cs
CA801FCH	250 lb.	5 ft. 0 in.–6 ft. 1 in. (Adult)	Forearm, push button, adjustable, aluminum	10 pr/cs

Crutch replacement parts

Cat. no.	Description	Packaging
CACP20028	Underarm pad, gray	100 ea/cs
CACP20028B	Underarm pad, heavy duty, gray	100 ea/cs
CAACK20026	Crutch tip, underarm pad and hand grip	50 pr/cs
CAACK20026B	Parts kit, heavy duty, crutch tip, underarm pad and hand grip	50 pr/cs
CAACT20027	Crutch tip, gray	100 ea/cs
CAACT20027B	Crutch tip, heavy duty, gray	100 ea/cs



CA901AD



CA801ADB



CA801FCH



CAACK20026B