



ComfortZone™ Single Socket Gel Liner

- 3 ply, 3mm of exposed gel, CoolMaxÆ exterior fabric, tapered shape
- Optimum cushion and protection
- Prevents abrasion and skin breakdown
- Unmatched patient comfort and satisfaction
- Not made with natural rubber latex & hypoallergenic

ITEM # - 12725 , 12755 , 12765 ,

1/Box - See next page for sizing details.



ComfortZone™ Double Socket Gel Liner

- 5 ply, 3mm gel
- Two layers of CoolMaxÆ fabric
- Tapered shape
- Prevents abrasion and skin breakdown
- Optimum cushion and protection
- Not made with natural rubber latex & hypoallergenic

ITEM # - 12625 , 12655 , 12665 ,

1/Box - See next page for sizing details.



ComfortZone™ Ultra Cushion Gel Liner





- Available in 3 ply or 5 ply, 3mm gel
- Two layers of cotton material
- Nylon-reinforced distal end
- Prevents abrasion and skin breakdown
- Optimum cushion and protection
- Not made with natural rubber latex & hypoallergenic

ITEM # - 12515 , 12525 , 12535 , 12555 , 12565 , 12575 ,

1/Box - See next page for sizing details.



COMFORTZONE™ SINGLE SOCKET GEL LINER SIZING GUIDELINES

ITEM#	SIZE	PLY	GEL LENGTH	LENGTH OF LINER	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
						
12725	Short/Narrow	3	10" (25cm)	12" (30.5cm)	9"-12" (23cm-30.5cm)	15"-18" (38cm-45cm)
12755	Short/Regular	3	10" (25cm)	12" (30.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)
12765	Medium/Regular	3	12" (30.5cm)	16" (40.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)

All sizes are approximate.

COMFORTZONE™ DOUBLE SOCKET GEL LINER SIZING GUIDELINES

ITEM#	SIZE	PLY	GEL LENGTH	LENGTH OF LINER	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
12625	Short/Narrow	5	10" (25cm)	12" (30.5cm)	9"-12" (23cm-30.5cm)	15"-18" (38cm-45cm)
12655	Short/Regular	5	10" (25cm)	12" (30.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)
12665	Medium/Regular	5	12" (30.5cm)	16" (40.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)

All sizes are approximate.

COMFORTZONE™ ULTRA CUSHION GEL LINER SIZING GUIDELINES

ITEM#	SIZE	PLY	GEL LENGTH	LENGTH OF LINER	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
12555	I-Short/Regular	3	10" (25cm)	12.5" (32cm)	7"-11" (18cm-28cm)	8"-13" (20cm-33cm)
12515	I-Short/Regular	5	10" (25cm)	12.5" (32cm)	7"-11" (18cm-28cm)	8"-13" (20cm-33cm)
12565	II-Medium/Regular	3	12" (30.5cm)	15"-17" (38cm-43cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
12525	II-Medium/Regular	5	12" (30.5cm)	17"-19" (43cm-48cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
12575	III-Medium/Wide	3	12" (30.5cm)	17"-19" (43cm-48cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)
12535	III-Medium/Wide	5	12" (30.5cm)	18"-21" (46cm-53cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)

All sizes are approximate.