



MATERNITY SIZE Pre-Pregnancy Dress Size	
Small	3 - 8
Medium	9 - 14
Large	15 - 18
X-Large	19 - 22
2X	23 - 26

Breathability	3
Compression	2
Flexibility	2
Stability	4

### Maternity Support

**Provides Extra Support and Helps Relieve Pressure on Abdomen**

Unlike other maternity supports, our Maternity Support creates an upward and inward lift, providing relief from strained muscles. By supporting the weight of the abdomen and transferring it to the spine, this support creates more comfort for both mom and baby. Elastic construction with soft, comfortable, hook and loop front panel to provide additional lift. Sized according to pre-pregnancy dress size. White. Available in S, M, L, XL, 2X, or 3X.

ITEM	DESCRIPTION
LSB-6090	Maternity Support



Wrap support around back and attach strap just below belly.

Attach front panel and stretch across lower abdomen. Attach to other side.

Readjust end, pulling upward and inward to lift abdomen.

### Better Binder™ Post-Partum Support Post-Operative, Post-Partum Back Stabilization

Designed by a physical therapist with over 40 years experience in women's health and orthopedics. The versatile and shape-adaptive Better Binder holds and hugs the post-operative belly, remodels lymphedemic bellies and gives bariatric patients comfort and control. Helps support back muscles and stabilizes excessive joint motion in the pelvis, spine or ribs. After pregnancy, Better Binder helps empower abdominal muscles and back, stabilizes excessive joint motion in the pelvis and remodels the belly.

Can be worn over or under clothing (but not directly over a wound). Latex-free elastic is soft but strong with perfect stretch. Edges won't cut in or fold over. Beige. Available in S, M, L or XL.



SIZE	WAIST
Small	20" - 29"
Medium	30" - 39"
Large	40" - 55"
X-Large	56" - 75"



Breathability	4
Compression	4
Flexibility	4
Stability	2

ITEM	DESCRIPTION
BBH-6906	Better Binder Post-Partum Support



MORE BABY HUGGER SELECTIONS AVAILABLE at [www.coreproducts.com/babyhugger](http://www.coreproducts.com/babyhugger)

MATERNITY SIZE Pre-Pregnancy Slack Size	
X-Small	0 - 4
Small	4 - 8
Medium	8 - 12
Large	12 - 16
X-Large	16 - 24
2X	24 - 28

*"The Baby Hugger is a godsend! I started wearing one when I was 32 weeks into my third pregnancy and never developed the nagging backaches I experienced with my first two."*

—Diana Brown, WI

### Baby Hugger® Maternity Support

Help relieve backache, pelvic pressure and heaviness, belly weakness and pain, hernias, swollen veins, hip pain, frequent urination or edema with our Baby Hugger Support. Our support helps protect ligaments from overstretching and causing prolapses and future incontinence. Baby Hugger's unique *central lift* takes weight off the back, bladder, veins and ligaments, and lifts thoroughly - from the top, beneath and at the sides. Cool, comfortable cotton fabric is completely adjustable, fitting for 9 months of pregnancy. Soft hook and loop closure. Crotch panel opens easily for bathroom use. Available in XS, S, M, L, XL or 2X.

ITEM	DESCRIPTION
BBH-6900	Baby Hugger Maternity Support Single
BBH-6902	Baby Hugger Extra Pant



BabyHugger®