

IV  
Solutions

Enteral  
Care

Respiratory  
Care

Medical/  
Surgical

Urology

The **Dale Bendable ArmBoard** positions and stabilizes the wrist to prevent unwanted movement when using arterial and peripheral lines in the critical care units, operating room, recovery room, and patient care units.

**Dale**<sup>®</sup>  
Always Better



## Bendable ArmBoard

THE ARMBEARD CAN BE CUSTOM-SHAPED TO ANY DESIRED POSITION.

The ArmBoard can be shaped to achieve more prominent radial artery access. A slight bend of the ArmBoard will provide the patient with a more natural hand/wrist position. The small and x-small sizes can be shaped to accommodate an infant's wrist or heel. The ArmBoards are covered in soft, comfortable material, and fastened to the wrist/heel with adjustable straps with hook and loop closure which allows the clinician easy access to the catheter site.



- Bendable inner core
- Comfortably maintains wrist in hyperextended position
- Secures the wrist
- Improves skin care
- Improves patient comfort
- Trouble-free monitoring
- Not made with natural rubber latex.



MR Conditional. See Conditions For Use by visiting our website for detailed information on scanning safely with this device.

### ORDERING INFORMATION

Product Number	Size	Quantity
650	Large, 23cm x 9cm (9" x 3.5")	10 per box
651	Medium, 13.7cm x 6.4cm (5.38" x 2.5")	10 per box
652	Small, 11.4cm x 4.4cm (4.5" x 1.75")	10 per box
653	X-Small, 10.8cm x 2.5cm (4.25" x 1")	10 per box

**Dale**<sup>®</sup>  
[www.dalemed.com](http://www.dalemed.com)

Dale<sup>®</sup> is a registered trademark of Dale Medical Products, Inc. in the USA & EU.  
©2017 Dale Medical Products, Inc.

PE-650-653 REV F