

DEBUG ALUMINUM ALL-TERRAIN WALKER



DEBUG ALUMINUM ALL-TERRAIN WALKER

debug

The Debug Aluminum All-Terrain Walker provides stability to individuals in walking on sandy beaches, rough terrain and even snow that normal walkers cannot accommodate. Light weight, portable and ready to go when you are.

LIGHT WEIGHT ALUMINUM FRAME FOLDS FOR TRANSPORTATION AND STORAGE / HAND BRAKE AND LOCKING MECHANISM / FOLDING SEAT AND BACK-PAD TO SIT / INCLUDED STORAGE BASKETS / STAINLESS STEEL CASTERS AND AXLES RESIST SALT AND SAND / PNEUMATIC WHEELS ROLL ACROSS SAND, SNOW AND GRAVEL / FITS THROUGH A STANDARD ADA DOORWAY / LOW MAINTENANCE DERLIN BRUSHINGS / ADJUSTABLE AND REMOVABLE PUSH HANDLES

SPECIFICATIONS

SPEC	ENGLISH (INCH)	METRIC (CM)
Overall height (adjustable)	34-39	86-142
Width	39	99
Weight	25 lbs	11 kg

ADDITIONAL NOTES

Available in weight capacities from 200 – 300 lbs

Frame material: anodized aluminum

Walker frames are available in pediatric, adult and large adult sizes

Casters and axles are made from #316L stainless steel with 10 year warranty

The Debug Rollator Style Walker is intended for Outside Use. It may not easily fit through doorways. While it is easily transportable, and folds, it is also heavier than standard rollators. Some users may require assistance with set up and transportation.

Manufactured in the USA, Debug has a 10 Year Warranty and a 30 day Money Back Guarantee.

© Copyright Deming Designs Inc., 2001–2016. All rights reserved.

ORDER ONLINE WWW.DEBUGPRODUCTS.COM