



FAWSSit™
PORTABLE SHOWERS

*Because people in wheelchairs
deserve warm, safe, easy to use,
accessible and affordable showers too!*

1-877-329-7748
www.fawssit.com

Would you like to improve the health and quality of life for a wheelchair-bound family member?

Don't want an expensive solution that permanently alters your home and requires permits, contractors and lots of time to implement?

FAWSsit™ showers come completely assembled and set up quickly in any room with water/drain access. Roll-in, shower and roll-out. Afterwards you simply fold the shower flat and stow away until the next use.

FAWSsit™ showers provide the full health benefits of a regular shower without need for lifts, or transfers.

Using the FAWSsit™ portable shower can decrease incidence of skin infections and UTI's, promotes circulation and muscle relaxation, and helps maintain critical body temperature control.

And showering IN the wheelchair is safe, without the dangers of slipping or dropping or falling.

The FAWSsit™ shower is lightweight and simple to use - no tools are required.

Makes any house, mobile home, or apartment completely wheelchair accessible for bathing. Great for travel too!



FAWSSIT SHOWERS COME IN FOUR SIZES TO ACCOMMODATE YOUR NEEDS



S2000 Standard (36" x 36" x 48")*

For wheelchair users who do not require a tilt back shower wheelchair. Carry case on wheels is also available.



SS4040 Super Standard (42" x 42" x 48")*

Same design as the S2000 but allows an extra 6" in width and depth for larger wheelchairs or footrests.

B5000 Bariatric (46" x 46" x 44")

For wheelchair users who require more space for larger size wheelchairs or special patient requirements.



R3000 Recliner (36" x 76" x 36")*

For the High quadriplegic, CP and any medical condition requiring a tilt back shower wheelchair.



* outside dimensions

BENEFITS ASSESSMENT GUIDE

Benefit	Conventional Bath				Shower		FAWSSit
	Chair	Side/Lift	Built-in	Walk in	Chair	Built in	
Active clean water flow					•	•	✿
Transfer from wheelchair not needed						•	✿
Ability to wash hair (without changing position)					•	•	✿
Ability to clean under skin folds					•	•	✿
Ability to soak feet	•	•	•	•			✿
No wait for tub to fill		•			•	•	✿
Perineal cleaning		•	•	•			✿
No splash over				•	•	•	✿
Caregiver stays dry	•	•	•	•			✿
Reclining position possible			•			?	✿
Easy/safe access to temperature controls				•		?	✿
Accessible to bariatrics						•	✿
Ability to “grow” with patient	•				•		✿
Ability to easily relocate and for travel	•				•		✿
COST	\$	\$\$-\$\$\$\$	\$\$\$	\$\$\$\$\$	\$	\$\$\$\$\$	\$\$\$

\$= Less than \$100 \$\$=Less than \$1000 \$\$\$= \$1000-\$3000 \$\$\$\$=\$3000-\$10,000 \$\$\$\$\$=\$10,000 - \$25,000 ? Depends on configuration

WHO BENEFITS FROM FAWSSIT SHOWERS:

ANYONE who has found traditional bathing difficult or unsafe.

SENIORS who have difficulty getting in and out of a shower or tub

PERMANENTLY DISABLED such as paraplegics, quadriplegics, or those with Multiple Sclerosis, Muscular Dystrophy, ALS, Spina Bifeda, Cerebral Palsy, Polio, or any condition which has dictated the need to use a wheelchair.

TEMPORARILY OR PARTIALLY DISABLED such as hip or knee replacement patients or partial amputees or the obese for whom standing is difficult or unsafe

CAREGIVERS who have had to risk their own health and comfort by lifting or bathing clients

FAMILY MEMBERS who are trying to avoid nursing homes for their loved ones

RENTERS OR FINANCIALLY CONSTRAINED INDIVIDUALS for whom re-modeling is not an option



FAWSsit™

PORTABLE SHOWERS

are extremely easy to setup without tools or any need for heavy lifting or straining



Setup for showering consists of unfolding the panels of the shower, and the bottom rubber pan, connecting the pump to electrical outlet and connecting the showerhead to faucet and the drain hose to sink drain.

Roll the chair with the patient and secure the bottom pan to the lowered support bar.



After showering, reverse the process to fold flat & store.



FAWSsit™ wheelchair showers are designed and built for long term use. From our customers we know that our earliest models are still working even after more than eight years of regular use.

They are 100% U.S. made, and feature a sturdy metal frame and custom designed, UL listed, mounted pump.

They come fully assembled, ready to use out of the box.

MedTrade® Atlanta Innovation Award Winner

FAWSsit WARRANTY INFORMATION:

10-year warranty on the anodized aluminum frame.
Full one year warranty on all other components.
Extended warranties are also available.

