

Leg Bag Systems



Hollister leg bags maintain a low profile as they rest discreetly against the leg with soft, pleated sides and are held comfortably in place with wide fabric straps.



Portrait:
Bastien Perret, France
"It's not possible for anyone to live life and never have any problems. But for every problem, there is a solution. A person just has to go out and find it."



Leg Bag Systems

The latex-free Hollister vented and non-vented leg bags use an anti-reflux valve to help prevent urine back flow. Our urinary leg bags are designed to be discreet and secure, promoting independence and quality of life.


Anti-Reflux Valve

- Helps prevent urine from back flowing into the tubing

Soft Vinyl Film

- Helps ensure comfort against the skin and limits unwanted odor

Optional Vented System



Vent

- Helps prevent air build-up in the leg bag, providing a lower profile.

Sterile




Convenience Drain

- Easy to operate with one hand, even for those with limited dexterity



Buttons

- Low profile buttons help secure the bag



Latex-free bedside drainage collection system with anti-reflux valve also available

Extension Tubing

- Our oval extension tubing helps keep the fluid pathway open
- Sterile and non-sterile options are available



Fabric Straps

- Wide straps help hold the bag comfortably and securely in place

Leg Bag



Urinary Leg Bag

Leg bag, latex-free

Stock No	Description	Length	Width	Box Qty
9814	Medium - 18 oz (540 mL)	10"	4-1/4"	10
9805	Large - 30 oz (900 mL)	11-1/8"	5"	10

Leg Bag Systems & Packs



Urinary Leg Bag Combination Pack

Includes:

- Leg bag
- Extension tubing (18 in) with connector, non-sterile
- 1 Pair leg bag straps – medium (23 in)

Latex-free; contains sterile components

Stock No	Description	Length	Width	Box Qty
9348	Medium (18 oz / 540 mL)	10"	4-1/4"	12
9349	Large (30 oz / 900 mL)	11-1/8"	5"	12



Vented Urinary Leg Bag Combination Pack

- Vented leg bag, sterile
- Extension tubing, kink-resistant, sterile
- 1 Pair leg bag straps – medium (23 in)

Latex-free; contains sterile components

Stock No	Description	Length	Width	Box Qty
9645	Medium (18 oz / 540 mL)	10"	4-1/4"	10
9655	Large (30 oz / 900 mL)	11-1/8"	5"	10

Accessories



Latex-Free Fabric Leg Bag Straps

Stock No	Description	Length	Width	Box Qty
9342	Small	15"	1-1/8"	10 Pair
9343	Medium	23"	1-1/8"	10 Pair



Latex-Free Extension Tubing with Connector

Stock No	Description	Length	Width	Box Qty
9345	Non-sterile	18"	–	10
9346	Sterile	18"	–	10

Bedside Drainage Collection System with Anti-Reflux Valve



Bedside Drainage Collection System with Anti-Reflux Valve

Sample Port
Sterile Fluid Pathway
Latex-Free

Stock No	Description	Length	Width	Box Qty
9839	68 oz (2000 mL)	12"	9-3/4"	20

A Hollister Technology

For detailed clinical questions concerning our products:
1.888.740.8999

For orders only:
1.800.323.4060

Hollister Incorporated
Libertyville, IL USA

Distributed by
Hollister Limited
Aurora, Ontario
1.800.263.7400

www.hollister.com



Hollister and logo, and Hollister are trademarks of Hollister Incorporated.

©2012 Hollister Incorporated. Printed in USA.
922093-1012