

Pressure-Activated Safety Lancet

McKesson Pressure-Activated Safety Lancets feature a one-step safety mechanism. Simply hold the lancet against the patient's fingertip and press down gently. The sterile needle will quickly penetrate skin to draw blood, then immediately retract back inside the safety mechanism. This helps ensure the entire needle is concealed before and after use in order to help minimize infection risk. These lancets are available in three gauge sizes with different penetration depths to provide optimal blood flow for different applications and skin types.

Features and Benefits

- Needle is concealed before and after use, minimizing exposure risk
- Advanced design helps mask the pain stimulus, resulting in greater comfort for the user while providing an accurate blood sample
- Innovative, single-use safety device helps prevent needle stick injuries and cross-contamination
- Sterile, tri-beveled needle for added protection
- 21G needle produces normal blood flow; ideal for standard capillary blood sampling
- 26G needle produces low blood flow; ideal for blood glucose and cholesterol tests, coagulation tests, and blood gas tests
- 30G needle produces micro blood flow; ideal for diagnostic testing that requires a small blood sample
- Not made with natural rubber latex



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Pressure-Activated Safety Lancet, 21g	16-PASL21G	100/bx, 20bx/cs
Pressure-Activated Safety Lancet, 26g	16-PASL26G	100/bx, 20bx/cs
Pressure-Activated Safety Lancet, 30g	16-PASL30G	100/bx, 20bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.