

Scrotal elastic truss
Model: 2936 (SINGLE), 2939 (DOUBLE)
HCPCS: L8300 (2936), L8310 (2939)



Effective in relieving discomfort while treating inguinal or scrotal hernias. This men's support has been designed to support most common reducible ruptures. The design of this truss allows the specially shaped scrotal pad to self-align and maintain its proper position, thus applying support as a constant, unified stable force.

- Measure around widest part of hips.
- SMALL 30.0-32.0 in.
 - MEDIUM 34.0-36.0 in.
 - LARGE 38.0-40.0 in.
 - X-LARGE 42.0-44.0 in.
 - 2X-LARGE 45.0-48.0 in.

For model 2936, select a garment for the side of the body affected (wearer's LEFT or RIGHT).

Lightweight single hernia support
Model: 2958
HCPCS: L8300



Very simple in design and very easy to use. Holds a single inguinal or scrotal hernia comfortably in reduction, and can also be used in special situations after a hernia repair to prevent recurrence. Its uncomplicated low profile asymmetric design allows it to be worn under clothing without detection. 4" web elastic encircling belt molds around body contours for a comfortable, custom fit.

- Measure around widest part of hips.
- SMALL 30.0-32.0 in.
 - MEDIUM 34.0-36.0 in.
 - LARGE 38.0-40.0 in.
 - X-LARGE 42.0-44.0 in.

Select a garment for the side of the body affected (wearer's LEFT or RIGHT).

Lightweight beige hernia support
Model: 2956
HCPCS: L8300, L8310



Plush lightweight material is contact closure compatible. Elastic waistband in back molds around body contours for a comfortable, custom fit. Easy bilateral hook and loop fastening and adjustment. Removable foam pads provide gentle, localized support at hernia site. Soft elastic understraps can be tightened or loosened to increase or lessen support. Completely free of metal buckles and snaps.

- Measure around widest part of hips.
- SMALL 30.0-35.0 in.
 - MEDIUM 35.0-40.0 in.
 - LARGE 41.0-46.0 in.

Lightweight bilateral hernia support
Model: 2958-D
HCPCS: L8310



Effective in relieving discomfort while treating inguinal or scrotal hernias. Often used after abdominal surgery to prevent or control incisional hernias. Adjusts easily to provide the ideal amount of compression. Removable foam pads provide localized support at hernia site. Detachable understraps can be tightened or loosened to increase or lessen support. Simple hook and loop fastening for easy application.

- Measure around widest part of hips.
- SMALL 30.0-32.0 in.
 - MEDIUM 34.0-36.0 in.
 - LARGE 38.0-40.0 in.
 - X-LARGE 42.0-44.0 in.