



DuraGel™ offers an unparalleled combination of durability, comfort, and skin protection for low to high activity level patients.

- Made with a durable gel formulation designed to help increase liner life and maintain comfort
- Ideal for moderately active patients engaging in every day living activities
- Nylon outer fabric engineered to provide superior abrasion resistance and gel bond
- Available in common pin and cushion profiles: 3mm, 6mm, and 9mm; Tapered and Uniform styles
- Packaged in a box



Pin Liner, 3mm Gel Thickness

Item #	Size	Style
18211	Medium	Uniform
18212	Medium +	Uniform
18213	Large	Uniform
18214	Large +	Uniform
18215	X-Large	Uniform

Item #	Size	Style
18216	Medium	Tapered
18217	Medium +	Tapered
18218	Large	Tapered

Pin Liner, 6mm Gel Thickness



Item #	Size	Style
18201	Medium	Uniform
18202	Medium +	Uniform
18203	Large	Uniform
18204	Large +	Uniform
18205	X-Large	Uniform

Item #	Size	Style
18206	Medium	Tapered
18207	Medium +	Tapered
18208	Large	Tapered

Pin Liner, 9mm Gel Thickness

Item #	Size	Style
18221	Medium	Uniform
18222	Medium +	Uniform
18223	Large	Uniform
18224	Large +	Uniform

DuraGel™ Sizing Guidelines

Proximal Circumference			Distal Circumference	
 Measure 15" from Distal End			 Measure 1.5" from Distal End	
Size	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" (23cm)	15" (38cm)	8" (20cm)	13" (33cm)
Medium+	13" (33cm)	20" (51cm)	9" (23cm)	13" (33cm)
Large	15" (38cm)	21" (53cm)	9" (23cm)	16" (41cm)
Large+	16" (41cm)	26" (66cm)	12" (30cm)	17" (43cm)
X-Large	16" (41cm)	29" (74cm)	12" (30cm)	21" (53cm)

All sizes are approximate.



Cushion Liner, 3mm Gel Thickness

Item #	Size	Style
18251	Medium	Uniform
18252	Medium +	Uniform
18253	Large	Uniform
18254	Large +	Uniform
18255	X-Large	Uniform

Item #	Size	Style
18256	Medium	Tapered
18257	Medium +	Tapered

Cushion Liner, 6mm Gel Thickness

Item #	Size	Style
18241	Medium	Uniform
18242	Medium +	Uniform
18243	Large	Uniform
18244	Large +	Uniform
18245	X-Large	Uniform

Item #	Size	Style
18246	Medium	Tapered
18247	Medium +	Tapered
18248	Large	Tapered

Cushion Liner, 9mm Gel Thickness

Item #	Size	Style
18261	Medium	Uniform
18262	Medium +	Uniform
18263	Large	Uniform
18264	Large +	Uniform