



SURFLO® Winged Infusion Sets

Our SURFLO Wing design minimizes the risk of vessel damage and enhances ease of penetration.

Compare our wings versus the competition and see why SURFLO is a highly respected brand.


Terumo's soft, flexible interlocking wings readily conform to the body's contour. The wing is angled for ease of positioning of the needle during venipuncture, for less risk of vessel damage or occlusion of the lumen.

The precisely honed needle with ultra-thin wall design assures ease of penetration and continued access with a single insertion.



These products contain no components made of natural rubber latex.

Surflo® Winged Infusion Sets



Product Code Description

3½" Tubing

Packed: 100/box, 5 boxes/case (500)

■ SV*19BLS	19G x ¾"	TW
■ SV*21BLS	21G x ¾"	UTW
■ SV*23BLS	23G x ¾"	UTW
■ SV*25BLS	25G x ¾"	TW

RW - Regular Wall
 TW - Thin Wall
 UTW - Ultra Thin Wall

Product Code Description

8" Tubing

Packed: 100/box, 10 boxes/case (1,000)

■ SV*25EL	25G x ½"	TW
■ SV*27EL	27G x ½"	RW

Product Code Description

12" Tubing

Packed: 50/box, 10 boxes/case (500)

■ SV-18BLK	18G x ¾"	TW
■ SV-19BLK	19G x ¾"	TW
■ SV-21BLK	21G x ¾"	UTW
■ SV-22BLK	22G x ¾"	UTW
■ SV-23BLK	23G x ¾"	UTW
■ SV-25BLK	25G x ¾"	TW

Call your Terumo representative at:

1.800.888.3786

or visit www.TerumoTMP.com

and click on Find Your Local Terumo Rep.

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURFLO is a trademark owned by Terumo Corporation, Elkton, MD. PM-00689



TERUMO MEDICAL PRODUCTS