

The New 3M™ Tegaderm™ Silicone Foam Border Dressing From the Innovative Leader in the Science of Securing to Skin



Gentle.
Absorbent.
Secure.

3M™ Tegaderm™ Silicone Foam Border Dressing easily handles high to low exuding wounds. It provides patients with comfort and clinicians with the confidence they're giving the best possible care.

Absorb & Evaporate

- Innovative fluid handling technology like no other on the market
- Unique four layer design provides high absorbency with breathability to reduce the risk of maceration
- Adapts to changing levels of exudate to maintain moisture balance for optimal wound healing

Comfort & Security


- 3M Silicone stays in place, yet removes gently making it ideal for fragile skin
- Soft, conformable and comfortable to wear
- Does not swell which decreases disruption of healing tissue
- Tegaderm™ film backing prevents strikethrough and protects against external bacteria and viruses*

*In vitro testing shows that the transparent film provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.



The New 3M™ Tegaderm™ Silicone Foam Border Dressing From the Innovative Leader in the Science of Securing to Skin


Ordering information:

Product Code	Overall Dressing Size	Foam Pad Size	HCPCS Code**
90640	3 in. x 3 in. 7.5 cm x 7.5 cm	1½ in. x 1½ in.	A6212
90641	4 in. x 4 in. 10 cm x 10 cm	2¾ in. x 2¾ in.	A6212
 90642	6 in. x 6 in. 15 cm x 15 cm	4¼ in. x 4¼ in.	A6213

**The reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coding for 3M products. Insurers' reimbursement policies can vary, however, and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level.

The Science Behind 3M's Excellent Fluid Handling Capacity

First, the soft foam layer rapidly wicks away exudate, which is then retained in the super absorbent layer. Next, the exudate is distributed across the moisture control layer before evaporating through the breathable film.



The 3M™ Tegaderm™ Foam Family of Products

3M™ Tegaderm™ Foam Dressings reduce the risk of maceration and desiccation by maintaining an optimal moisture balance so you can be sure of providing your patients with the best possible care. Select Tegaderm™ Silicone Foam Border Dressing for patients with fragile skin. Consider Tegaderm™ Adhesive Foam Dressing for the added benefit of wear times up to seven days. Both dressings provide excellent fluid handling.

For more information about the 3M Foam family of products, visit www.3m.com/Tegaderm-Foam, contact your 3M Critical and Chronic Care Solutions representative or call the 3M Health Care Customer Helpline at 1-800-228-3957. Outside of the United States, contact the local 3M subsidiary.

ValidityMedical.com
Add Vitality to Your Life

[DISCOUNT Tegaderm Dressings ON SALE at Validity Medical](#)



3M Critical & Chronic Care Solutions Division
2510 Conway Ave.
St. Paul, MN 55144
USA
1-800-228-3957
www.3m.com/C3SD

3M and Tegaderm are trademarks of 3M.
Please recycle. Printed in U.S.A.
© 3M 2013. All rights reserved.
70-2010-9225-4

Tegaderm™ Foam Silicone
Part of the Full Line of Tegaderm™ Brand Advanced Wound Care

Tegaderm™ High Performance
Foam Adhesive Dressing