

# Medical Leg Wear Pantyhose

## Standard Figure

● 1756 20-30 mmHg


## Maternity Figure

● 1757 20-30 mmHg



## Full Figure

● 1758 20-30 mmHg



## Material Information

1756 | 1757 | 1758

Manufactured using 80% nylon and 20% spandex

## Product Description

Designed to help relieve moderate conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. Color options are basic neutrals for all-occasion wear. For men and women.

## Product Features

- Therapeutic graduated compression
- Firm support.
- Classic, concealing smooth-knit texture fabric
- Flat panty seams for a comfortable inconspicuous fit
- Ventilated cotton panel for comfort
- Available with maternity panel (1757)

LATEX FREE

## Available Colors





\*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

## Classic Medical Compression Pantyhose

Model	Compression	Description	Latex Free
1756	20 - 30 mmHg	Standard Figure	✓
1757	20 - 30 mmHg	Maternity Figure	✓
1758	20 - 30 mmHg	Full Figure	✓



# Anti-Embolism

