

#8550/#8550-SH Recovery Couches



ISO 9001-2000 CERTIFIED

Perfect Low-Height Table for Rest and Recovery

Winco's recovery couches are a great choice for almost any medical facility or school nurse's office. With a height of only 19" and a weight capacity of 400lbs. a patient can sit or lie comfortably.



#8550
Recovery Couch



#8550-SH
Recovery Couch
w/ Shelf

#8550/#8550-SH Recovery Couches




ISO 9001-2000 CERTIFIED

Comfortable Medium Density Foam for Superior Comfort

Heavy Gauge Expanded Vinyl

High Strength Powder-Coated Steel Frame with a 400 lb. Weight Capacity



Adjustable Leveling Feet

Model #8550-SH Shown

Standard Paper Holder (Paper Cutter Not Shown)

VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Black	Burgundy	Mauve	Blueridge
02	03	04	06	07	08	14	16	17
	32	33	37	60	61	62	63	
	Grotto Optional	Eucalyptus Optional	Navy Optional	Pyrite Optional	Gypsum Optional	Ice Blue Optional	Cypress Optional	

SPECIFICATIONS #8550

OVERALL HEIGHT:	19" (48cm)
OVERALL WIDTH:	28" (71cm)
OVERALL LENGTH:	72" (183cm)
WEIGHT:	73lbs. (33.1 kg)*
WEIGHT CAPACITY:	400lbs. (181.4 kg)
FOOTPRINT:	28" x 72" (71cm x 183cm)

* Shelf 25lbs. Shelf Option

OPTIONS

Table Options

No Option Chosen = Plain Table
SH-Shelf

Height Options

N/A