


CLOG ZAPPER* ENTERAL FEEDING TUBE DECLOGGING SYSTEM

IN CLINICAL USAGE CLOG ZAPPER WAS SUCCESSFUL IN RESTORING PATENCY TO FEEDING TUBES CLOGGED WITH ENTERAL FORMULAS OVER 75% OF THE TIME ^{6 7}

Kit Includes:

- 1 x 10ml Oral syringe with declogging powder
- 1 x 10ml oral syringe (not used with ENFit® connectors on feeding tubes)
- 1 x 31cm applicator
- 2 x 10cm gauze pads.
- ***ENFit® syringe will be required for additional feeding tube water flushes***

- CLOG ZAPPER* is food graded and combines an “enzyme cocktail,” acids, buffers, anti-agents, and metal inhibitors.
- Ingredients all work together to loosen, break down and dislodge clogs.
- Enzymatic formula reduces clogging more effectively than water flushing alone.



CLOG ZAPPER* - Feeding Tube Declogging Kit

| Product Code | Product Description | Packaging |
|--------------|---------------------|-----------|
| 20-0002 | Clog Zapper* | 10/Case |

References:

- Frankel EH, Enow NB, Jackson KC, Kloiber LL. Methods of restoring patency to occluded feeding tubes. Nutr Clin Pract. 1998;13(3):129-131
- Lord LM. Restoring and Maintaining Patency of Enteral Feeding Tubes. Nutr Clin Pract. 2003;18(5):422-426