

Invacare®  
Rollators  
& Dolomite  
Rollators

Invacare® Rollators & Dolomite® Rollators



**E**xpand your mobility while adding comfort to your daily activities with an Invacare Rollator or a Dolomite Rollator.

Invacare Rollators include all the features you need to make walking easy and safe, like simple-to-operate, ergonomic hand brakes; built-in handles for easy folding; and a locking mechanism that prevents accidental folding. Our four-wheel models also offer adjustable handle height to ensure a custom fit, plus removable baskets and trays.

Dolomite Rollators can improve the mobility of people with many different conditions, reducing their risk of falling and increasing their independence and self-esteem. A wide range of rollators are available combined with unique accessory offerings.

If you are living with any of the following conditions, Invacare recommends the appropriate rollator.

#### **FRACTURED HIP OR HIP REPLACEMENT**

Recommended Rollator:

- All models, depending on client's preference and size

#### **STROKE PATIENTS**

Recommended Rollator:

- Alpha with one-hand brake, Symphony, with one-hand brake and HP handle

#### **RESPIRATORY CONDITIONS**

Recommended Rollator:

- All models (except Melody, Soprano) fitted with the Oxygen tank holder

#### **ARTHRITIS**

Recommended Rollator:

- Alpha Basic/Alpha Advanced with Forearm Supports, & models that can be fitted with Anatomical Handles: Symphony, Legacy, Maxi, Amigo, Spartan, Sandstorm, Frontier Lite, Stingray

#### **OSTEOPOROSIS**



## Dolomite Soprano Rollator

Model no. I2030-41-21 - Colour: Blue

The new Dolomite Soprano is a heavy-duty walker with large 10" wheels, providing excellent performance on rough, uneven ground, making it ideal for active outdoor users. With the unique cable-free braking system the Dolomite Soprano has a wide, sturdy frame and reinforced seat, allowing for an increased weight capacity of 355 lbs.

### FEATURES

- Cable-Free Braking System
- Frame – Sturdy heavy-duty frame provides excellent stability
- Adjustable Handles – Infinite adjustments to the required height and adjustable inward angle up to a maximum of 15°
- Hand Grips – Gently wedge-shaped hand grips provide a comfortable hand position for all hand sizes
- Flip-up Seat – Extra large seat flips up for greater reach access
- Folding – Stands when folded for easy storage and transport
- 10" wheels – large solid wheels provide excellent friction and a comfortable ride
- Basket and back support included.

### Specifications

	Soprano
Distance between rear wheels	22"
Minimum handle height	31"
Maximum handle height	39"
Seat Height	24"
Overall width	25 1/2"
Overall height	38"



## Dolomite Melody Rollator

Model no. I2010-42-21 MELODY Standard, Seat Height 24"

Model no. I2012-42-21 MELODY Low, Seat Height 21"

Model no. I2014-42-21 MELODY Super Low, Seat Height 19"

Colour: Cranberry Red

The Dolomite Melody is ahead of its time with its revolutionary braking system and advanced features. Safe, reliable and very easy to operate, the Dolomite Melody can be adjusted in many ways, making it suitable for most clients' needs.

### FEATURES

- Cable-Free Braking System
- Adjustable Handles – Infinite adjustments to the required height and adjustable inward angle up to a maximum of 15°
- Hand Grips – Gently wedge-shaped hand grips provide a comfortable hand position for all hand sizes
- Flip-up Seat – Extra-large seat flips up for easier reach access and gait training
- Folding – Stands when folded for easy storage and transport
- 8" Wheels – Long-lasting SLII wheels provide excellent friction, stainless steel brake heel prevents dirt build-up
- Basket and back support included.

### Specifications

	Super Low	Low	Standard
Distance between rear wheels	19"	21"	21"
Minimum handle height	26"	28"	31"
Maximum handle height	31"	35"	39"
Seat Height	19"	21"	24"
Overall width	22"	24"	24"
Overall height	25"	26"	27"



## Dolomite Symphony Rollator

Model no. I2160 SYMPHONY Standard  
 Model no. I2161 SYMPHONY Low  
 Model no. I2162 SYMPHONY Super Low  
 Colour: Silver

Made of ultra-light aluminum, the Dolomite Symphony represents the ultimate in safety, comfort and ergonomic design. Finished in an elegant silver colour, the Dolomite Symphony combines style and sophistication with the weight advantage of aluminum.

### FEATURES

- Lightweight aluminum frame
- Dolomite Advanced Brake System (DABS II)
- Flip-up seat for greater reach access and gait training
- Adjustable back support provides comfort and safety
- Infinite handle height and angle adjustments
- Back support and basket included
- Anatomical handles help prevent pain and numbness in the hands.

## Specifications

	Super Low	Low	Standard
Distance between rear wheels	18 1/2"	18 1/2"	18 1/2"
Minimum handle height	25"	27 3/4"	30"
Maximum handle height	30"	33 3/4"	36"
Min. & Max. tall handle heights	N/A	N/A	N/A
Seat Height	18"	21 1/2"	24 1/2"
Overall height	22"	22"	22"



## Dolomite Legacy Rollator

Model no. I2050-xx-xx LEGACY Standard

Model no. I2052-xx-xx LEGACY Low

Model no. I2054-xx-xx LEGACY Super Low

Colour: Candy Apple Red, Navy Blue

The Dolomite Legacy is a high quality indoor/outdoor rollator. Reliable, durable and user-friendly with more features and more value.

### FEATURES

- Flip-up seat for greater reach access and gait training
- Unique frame design with handles facing forward to promote better posture
- Patented brake system with adjustable tension
- Large 8" SLII wheels with superior bearings for increased maneuverability
- Tall handle model available
- Basket and adjustable back support included.

### Specifications

	Super Low	Low	Standard
Distance between rear wheels	19"	19"	19"
Minimum handle height	28"	31"	33"
Maximum handle height	31"	34"	37"
Min. & Max. tall handle heights	N/A	34-37"	35-39"
Seat Height	19"	21 1/2"	24"
Overall Height	22"	22"	22"



## Dolomite Amigo Rollator

Model no. I2360-36-2586 AMIGO Standard, Seat Height 24 1/2"

Model no. I2350-36-2586 AMIGO Low, Seat Height 21 1/4"

Colour: Hunter Green

The Dolomite Amigo is a sturdy, fully featured rollator suitable for both indoor and outdoor use. Designed for safety and comfort, the Dolomite Amigo offers an affordable answer to independent mobility and an active lifestyle.

### FEATURES

- Seat flips up for greater reach access and gait training
- Height adjustable handles
- Patented brake system with adjustable tension
- Basket and adjustable back support included.

## Specifications

	Low	Standard
Distance between rear wheels	16 3/4"	16 3/4"
Minimum handle height	27 1/2"	30 3/4"
Maximum handle height	34"	37"
Seat Height	21 1/4"	24 1/2"
Overall width	23 1/2"	23 1/2"
Overall height	27 1/2"	27 1/2"

# Dolomite® Rollators



## Dolomite Alpha Advanced/Basic Rollator

Model no. I2220-xx-xx ALPHA Advanced

Model no. I2230-xx-xx ALPHA Basic

Colour: Sapphire Blue

The Dolomite Alpha is a special rollator designed specifically for users that have to deal with the added complexity of walking with arthritic or stroke conditions. The unique, adjustable forearm supports and the braking system allows for maximum stability with minimum braking effort.

### FEATURES

- Frame construction maximizes stability and encourages proper posture and balance
- Forearm supports adjustments in angle, width and depth
- Individually adaptable handgrip can be adjusted in height and laterally in angle
- Unique brake grip is extremely easy to use and allows for limited hand dexterity
- Flip-up seat for greater reach access and gait training
- Dolomite Advanced Brake System (DABS II)
- Basket included.

## Specifications

	Advanced/Basic
Distance between rear wheels	24 1/2"
Minimum handle height	37 1/2"
Maximum handle height	46 3/4"
Seat Height	25 1/2"
Overall width	27 1/2"
Overall depth	27 1/2"





## Dolomite Maxi Rollator

Model no. I2190-xx-xx MAXI Standard

Model no. I2191-xx-xx MAXI Low

Colour: Sea Foam Green, Blue

The Dolomite Maxi is an extra wide rollator with a 400 lb. weight capacity and 20 1/2" of width between the handles, making it ideal for the heavier user. Choose the Dolomite Maxi for strength, durability, quality and style. Available in standard and low heights.

### FEATURES

- Reinforced frame with 400 lb. weight capacity
- Extra wide between handles
- Infinite handle height and angle adjustments
- Anatomical handles help prevent pain and numbness in the hands
- Dolomite Advanced Brake System (DABS II)
- Basket and adjustable back support included
- Excellent for both indoor/outdoor use.

### Specifications

	Low	Standard
Distance between rear wheels	21 1/2"	21 1/2"
Minimum handle height	28"	32"
Maximum handle height	34"	38"
Min. & Max. tall handle heights	N/A	N/A
Seat Height	21 1/2"	25 1/2"
Overall Width	27"	27"

## Dolomite® Rollator Accessories

An extensive selection of unique, functional accessories to meet many different needs.



**PADDED BACK SUPPORT:**

Made of soft velvet, for more comfort when seated.

Model no. I2514



**PADDED SEAT CUSHION:**

Foam cushion with washable cover is easy to attach to the rollator seat.

Model no. I2518

**CANE HOLDER:**

Easy to install and operate, is suitable for any single point cane and can be mounted on the front or rear of any of the rollators.

Model no. I2606



**OXYGEN TANK HOLDER:**

To meet respiratory equipment needs, this durable holder is suitable for "D" aluminum tank.

(Not available on Super Low)

Model no. I2682-25



**BAG FOR BASKET:**

The water-resistant, zippered bag attaches easily to the basket – extra security for personal belongings.

Model no. I2513



**FOOD AND BEVERAGE TRAY:**

Ideal for in-home use, this affordable and easy-to-clean tray simply attaches to the rollator seat.

Model no. 7031 I

**ALPHA PLATFORM:**



### **SLOW DOWN BRAKE:**

Unique design allows user to control speed for increased safety and gait training without wearing on the wheel. Ideal for Parkinson's or where cognition is an issue.

Model no. I2535-21 or I3041-21



### **(HP) HEMIPLEGIC HANDLE:**

The foam handle can be angled to suit individual needs to promote a safe grip on the weak side. Brake handle can be rotated for comfort. Can be mounted on right or left.

Model no. I2611-L-60, I2611-R-60



### **ANATOMICAL HANDLES:**

The larger handle surface reduces fatigue through improved weight distribution and prevents pain and numbness in the hands. Ideal for users with arthritis.

Model no. I2509



### **ONE-HAND BRAKE:**

If hand function is impaired on one side, the one hand brake allows for braking or locking of both brakes from the user's strong side. Ideal for arthritic or stroke clients.

Model no. I2505-L, I2505-R



### **TALL HANDLES:**

These chrome handles with durable foam grips improve rollator fit and comfort giving extra handle height for those who need it.

Model no. I3235-AK

## **Accessories Availability**

	Soprano	Melody	Symphony	Legacy	Maxi	Alpha	Amigo
Slow Down Brake			•	•	•	•	
One-Hand Brake			•	•	•	•	•
Tall Handles				*•			•
Anatomical Handles			•	•	•		•
HP Handle			•				

## Invacare® Rollite™ Rollators

Model no. 65100 - Rollite Rollator w/5" wheels Adult - Colour: Midnight Blue

Model no. 65100-JR - Rollite Rollator w/5" wheels Junior - Colour: Midnight Blue

Model no. 65100R - Rollite Rollator w/5" wheels Adult - Colour: Electric Red

Model no. 65100R-JR - Rollite Rollator w/5" wheels Junior - Colour: Electric Red

Model no. 68100 - Rollite Rollator w/8" wheels Adult - Colour: Midnight Blue

Model no. 68100-TA - Rollite Rollator w/8" wheels Tall Adult - Colour: Midnight Blue

The Invacare line of Rollite™ Rollators is unlike any other in the industry. Invacare has taken the Rollite one step further, offering a new line which addresses the sizes, wheel types, and colours that you've been looking for. The new Rollite line offers adult and junior sizes in different colours. For improved mobility over rougher terrain, large 8" wheel models are available in adult and tall sizes (one colour). Along with additional features, the new Rollite line provides the comfort, durability and quality that you've come to expect from Invacare products.

### FEATURES

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height.



## Specifications

	65100 Rollite Adult (Midnight Blue)	65100-JR Rollite Junior (Midnight Blue)
Depth opened at base	24"	24"
Depth folded	7"	7"
Width inside handgrips	18"	18"
Width inside base of legs	22"	22"
Width overall	27"	27"
Height adjustment min.-max.	32-36"	29-33"
*Approx. client height min.-max.	5'4"-6'	4'11"-5'7"
Seat-to-floor height min.-max.	19-23"	16 1/2"-20 1/2"