

# The Leader in Compression Therapy



**Post-operative  
Compression  
Garments**

#1 physician recommended

**JOBST**®

*Comfort, Health and Style!*

## Plastic Surgery Girdle

- Delivers effective compression following liposuction to reduce swelling and ensure optimal results
- Provides support after abdominal, buttock, flank or thigh surgery
- Latex free
- Made with open-weave, powernet fabric to help keep patients cool and dry
- Non-binding elastic holds and helps prevent rolling
- Open pubis
- Waist measurement determines size
- Mid-thigh and long-leg styles available
- Large break-resistant zippers make it easy to put on unassisted



Open pubis front



Open pubis back

## Epstein™ Facioplasty Support & JOBST® Faci



Epstein™ Facioplasty Support



JOBST® Facioplasty Elastic Support

## Surgical Vest

The JOBST® Surgical Vest provides increased patient comfort and satisfaction. Patient comfort is a concern during post-operative recovery, as it may affect wearing compliance and recovery time.

- Features break-resistant nylon zipper and VELCRO® hook and eye closures in the front for easy adjustments
- Supports breast and chest tissue after breast surgery
- Helps manage post-operative edema
- Luxuriously soft fabric
- Latex free
- Minimizes hematoma formation
- Holds surgical implants in place
- Machine washable and dryable
- Available with or without cups
- Reduces the need for surgical dressings
- Allows easy examination of the surgical site
- Helps to reduce tension on surgical incision



Front fastening Velcro® hook and eye closures



Break-resistant nylon zipper



New luxuriously soft fabric



## plasty Elastic Support

- Provides effective pressure and support following facial surgery and face lifts
- Aids in controlling post-operative edema
- Holds dressings or moist compresses in place without tape
- Adjustable Velcro® brand closures
- Easy for patient to put on without assistance
- JOBST® Facioplasty can be slit to expose patient's ears
- JOBST® Facioplasty is made of a latex free fabric that is comprised of a blend of nylon and spandex
- Epstein Facioplasty is made with beige tricot fabric

JOBST® Surgical Vest without cups



# Ordering Information

## JOBST® PLASTIC SURGERY GIRDLE

### Mid-Thigh (Female)

Fits Hips 8" - 11" (20-28 cm) Larger Than Waist

110661	Size Small, Waist measurement 24" - 26" (61-66 cm)
110662	Size Medium, Waist measurement 27" - 28" (69-71 cm)
110663	Size Large, Waist measurement 29" - 30" (74-76 cm)
110664	Size Extra Large, Waist measurement 31" - 32" (79-81 cm)
110670	Size 1 Extra Large, Waist measurement 33" - 34" (84-86 cm)
110672	Size 2 Extra Large, Waist measurement 35" - 36" (89-91 cm)

### Long-Leg (Female)

Fits Hips 8" - 11" (20-28 cm) Larger Than Waist

110666	Size Small, Waist measurement 24" - 26" (61-66 cm)
110667	Size Medium, Waist measurement 27" - 28" (69-71 cm)
110668	Size Large, Waist measurement 29" - 30" (74-76 cm)
110669	Size Extra Large, Waist measurement 31" - 32" (79-81 cm)
110671	Size 1 Extra Large, Waist measurement 33" - 34" (84-86 cm)
110673	Size 2 Extra Large, Waist measurement 35" - 36" (89-91 cm)

### (Male)

Fits Hips 4" - 7" (10-18 cm) Larger Than Waist

110680	Size Small, Waist measurement 30" - 32" (76-81 cm)
110681	Size Medium, Waist measurement 34" - 36" (86-91 cm)
110682	Size Large, Waist measurement 38" - 40" (96-102 cm)
110683	Size Extra Large, Waist measurement 42" - 44" (107-112 cm)
110684	Size 1 Extra Large, Waist measurement 46" - 48" (117-122 cm)
110685	Size 2 Extra Large, Waist measurement 50" - 52" (127-132 cm)

## EPSTEIN™

### Facioplasty Support for Neck and Chin

110542	One Size Fits All
--------	-------------------

Please refer to the product label and/or package insert for full instructions on the safe use of these products.  
VELCRO® is a registered trademark of VELCRO USA, INC.

## JOBST® SURGICAL VEST

### With Cups Sizing

Size 1: 30C, 32B, 32C, 34A, 34B, 36A  
Size 2: 34C, 34D, 36B, 36C, 38A, 38B  
Size 3: 34DD, 36D, 36DD, 38C, 38D, 40C  
Size 4: 38DD, 40D, 40DD, 42C, 42D, 44C  
Size 5: 42DD, 44D, 44DD, 46C, 46D, 48C

### Without Cups Sizing

Chest circumference (under arm measurement)

Size 1: 31" - 35" (79-89 cm)  
Size 2: 35 1/8" - 39" (89-99 cm)  
Size 3: 39 1/8" - 43" (99-109 cm)  
Size 4: 43 1/8" - 47" (109-119 cm)  
Size 5: 47 1/8" - 51" (120-130 cm)

### With Cups

111901	Size 1
111902	Size 2
111903	Size 3
111904	Size 4
111905	Size 5

### With Left Cup Only

111911	Size 1
111912	Size 2
111913	Size 3
111914	Size 4
111915	Size 5

### Without Cups

111906	Size 1
111907	Size 2
111908	Size 3
111909	Size 4
111910	Size 5

### With Right Cup Only

111916	Size 1
111917	Size 2
111918	Size 3
111919	Size 4
111920	Size 5

## JOBST® Surgical Vest

with cups

Size 1	30				C		
Size 2	32		B		C		
Size 3	34	A	B	C	D	DD	
Size 4	36	A	B	C	D	DD	
Size 5	38	A	B	C	D	DD	
	40			C	D	DD	
	42			C	D	DD	
	44			C	D	DD	
	46			C	D		
	48			C			

## FACIOPLASTY

### Elastic Support for Ears, Cheek and Chin

Measurements determined by circumference from tip of chin to crown of head

111825	Size Small, under 25" (64 cm)
111826	Size Medium, 25" - 26 7/8" (64-68 cm)
111827	Size Large, 27" or more (69 cm)



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633  
Tel. (+1) 704 331 0600 Fax (+1) 704 358 4558 www.bsnmedical.com  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.jobst.com  
To order toll-free: BSN medical (+1) 800 552 1157; JOBST (+1) 800 537 1063

JOBST - a brand of BSN medical



Comfort, Health and Style!



## JOBST® Donning Aids and Accessories



• Easy to Use!

Donning Devices help the patient put on JOBST® compression stockings  
Accessories are handy tools that provide convenience and value



Comfort, Health and Style!