

Cleansing/Irrigation

Sitz Baths

- Fitting standard toilets and commodes, sitz baths are portable and can be used in a variety of healthcare settings (Hospital, Long Term Care, and Home Care)
- With rounded edges for comfort and ease of handling, sitz baths are perfect for administering soothing irrigation directly to the perineal area
- Each set contains a 200 mL bag, 60" tubing, adjustable clamp, basin, and easy to follow instructions printed directly on the bag
- H995 elongated sitz bath is designed for use with water saving, low water consumption (6.0 Lpf/1.6Gpf) style toilets with elongated bowl design



H990

Code#	Color	L x W x FH x BH (inches)	Qty/cs	Packaging
2000 mL Bag w/60" Tubing, (Material: Polypropylene Basin and Clamp / PVC Bag and Tubing).				
H990-05	Gold	15 x 15 x 4.7 x 4.3	10	Individually Bagged
H990-07	Turquoise	15 x 15 x 4.7 x 4.3	10	Individually Bagged
H990-10	Dusty Rose	15 x 15 x 4.7 x 4.3	10	Individually Bagged
Elongated, 2000 mL Bag w/60" Tubing, Graduated: 500 mL increments up to 2000 mL and 1 and 2 qt measurements.				
H995-10	Dusty Rose	16 x 13 x 4.6 x 4.9	14	Individually Bagged



H995

Cleansing Enema Set

- Convenient for administering enemas, cleansing sets place necessary supplies within easy reach
- Choose from two types of enema administration sets — bucket or bags
- Both hold 1500 mL of fluid and include 60" tubing with pre-lubricated tip, adjustable clamp, castile soap, waterproof drape, and easy to read instructions

Code#	Capacity	Qty/cs	Qty/pkg
Graduated: 500 mL increments up to 1500 mL, (Material: Polypropylene Bucket / PVC Tubing).			
2560	1500 mL Bucket/60" Tubing	50	Polybag
Graduated: 100 mL and 4 oz increments up to 1500 mL or 48 oz, (Material: PVC Bag and Tubing).			
2562	1500 mL Bag/60" Tubing	50	Box



2560



2562

The material used for all the above plastic products is polypropylene unless otherwise indicated.