



## Shower Chair with Back and Removable Padded Arms

- Removable padded armrests provide comfort and added security
- Blow-molded bench and back provide comfort and strength
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Leg height adjusts in 1" increments

PRODUCTS	Description
12445KD-1	With Knock Down Back, 1/cs
12445KD-2	With Knock Down Back, 2/cs

SPECIFICATIONS	Description
Outside Legs	17.75 (W) x 19.25 (D)
Seat Width	16"
Seat Depth	14.5", 16.25"
Seat Height	16" - 21"
Weight (-1   -2)	9.75 lbs.   19 lbs.
Carton Shipping Weight (-1   -2)	10 lbs.   21 lbs.
Width Between Arms	18.5"
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
12002T	Replacement Tip 4/bx