


Multi-Ply ShearCare1500 Pressure Redistribution Foam Mattress



1500SC-1-FB

- Bariatric, latex free, therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Full width top layer is high-resilient, high-density die-cut foam for pressure redistribution and maximum support
- High-density foam bottom layer provides support and stability for the bariatric patient
- Patented LuraQuilt foam heel section offloads pressure from vulnerable heels
- Removable premium, multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS

PRODUCTS	Description
1500SC-42-FB	42" Multi-layered/Multi-zoned Bariatric Foam Mattress, E0184, 1/ea
1500SC-48-FB	48" Multi-layered/Multi-zoned Bariatric Foam Mattress, E0184, 1/ea

SPECIFICATIONS

SPECIFICATIONS	Description
Cover Material	Fluid Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	2
Side Rails	No
Dedicated Heel Section	Yes
Weight Capacity	42" 650 lbs., 48" 750 lbs.
Warranty	5 Year Limited Foam 1 Year Cover

ACCESSORIES

ACCESSORIES	Description
CPU64280	Replacement Cover, 1/ea
CPU64880	Replacement Cover, 1/ea
CPU65480	Replacement Cover, 1/ea