


### Universal Suction Tubing and Filter Kit



HCPCS  
A7002

Item #	Product Description	UOM
18600-KITN	Universal Suction Tubing and Filter Kit	1/bx

- Includes hydrophobic suction filter, 10" soft flexible medical grade silicon suction tubing, 72" suction tubing with blue tip, universal suction canister elbow.

### Hydrophobic Filter for Suction Machines




Item #	Product Description	UOM
18600-FILTER	Hydrophobic Filter for Suction Machines	3/bx

- Replacement hydrophobic filter compatible with most suction machines and tubing on the market.

### Suction Tubing



SUCP TUBING 10



SUCP TUBING 72



HCPCS  
A7002

Item #	Product Description	UOM
SUCP TUBING 10	Suction Tubing Blue Tip 10"	1/ea
SUCP TUBING 72	Suction Tubing Blue Tip 72"	1/ea

**COMPARE**  
**SUCP TUBING 72 TO**  
DeVilbiss® 7305D-612

### Universal Flexible Elbow for Suction Machines



TO CANISTER



Item #	Product Description	UOM
18600-ELBOW	Universal Flexible Elbow for Suction Machines	1/ea

- Replacement universal flexible elbow.
- Connects suction canister directly to filter, eliminating the need for tubing connections.

**COMPARE TO**  
DeVilbiss® 73050-619