


CARESITE® (LAD) System

Description	Product Code	REF Number	Unit Per cs
-------------	--------------	------------	-------------

CARESITE Standard Bore Extension Sets

Components on this page are not made with DEHP or natural rubber latex.




Standard bore extension set with non-removable slide clamp, CARESITE valve injection site, SPIN-LOCK® connector. Priming volume: 1.5 mL, Length: 8 in. (20 cm)

354219

354219

50




Standard bore extension set with two CARESITE valve injection sites 28 in. and 6 in. above distal end, non-removable slide clamp, SPIN-LOCK connector. Priming volume: 5.5 mL, Length: 35 in. (89 cm)

354220

354220

50

Filterflow® Filtered Extension Sets with CARESITE Injection Site




0.2 micron filter, CARESITE valve injection site 6 in. above distal end, on/off clamp, SPIN-LOCK connector. Priming volume: 5 mL, Length: 16 in. (40.6 cm)

354221

354221

50




1.2 micron filter, CARESITE valve injection site 6 in. above distal end, on/off clamp, SPIN-LOCK connector. Priming volume: 5 mL, Length: 16 in. (40.6 cm)

470117

470117

50

Filterflow Filtered Administration Sets with CARESITE Injection Sites



(15 drops/mL) Universal spike, backcheck valve, CARESITE injection sites 87 in. and 6 in. above distal end, roller clamp, 0.2 micron filter, slide clamp, SPIN-LOCK connector. Priming volume: 19 mL, Length: 104 in. (264.2 cm)

354207

354207

50



(15 drops/mL) Universal spike, roller clamp, 1.2 micron filter, slide clamp, CARESITE injection site 6 in. above distal end, SPIN-LOCK connector. Priming volume: 23 mL, Length 98 in. (248.9 cm)

470115

470115

50