







ULTRASITE® Valve Needleless System

Description	Product Code	REF Number	Unit Per cs
Standard Bore Extension Sets with ULTRASITE Injection Site			
 <p>Standard bore extension set with ULTRASITE valve, non-removable slide clamp, SPIN-LOCK® connector. Components are not made with DEHP or natural rubber latex. Priming volume: 0.7 mL, Length: 6.5 in. (17 cm)</p>	CSE8SL	473444	100
 <p>Standard bore extension set with bonded ULTRASITE valve, non-removable slide clamp, SPIN-LOCK connector. Components are not made with DEHP or natural rubber latex. Priming volume: 1.0 mL, Length: 6 in. (15 cm)</p>	CSE8SLB	473447	100
 <p>Extension set with bonded ULTRASITE valve, kink-resistant tubing, non-removable slide clamp, SPIN-LOCK connector. Components are not made with natural rubber latex. Priming volume: 0.68 mL, Length: 7 in. (18 cm)</p>	US1320	474921	50
 <p>Standard bore T-port extension set with ULTRASITE valve, non-removable slide clamp, SPIN-LOCK connector. Components are not made with DEHP or natural rubber latex. Priming volume: 1.3 mL, Length: 7.5 in. (18.8 cm)</p>	CSE8TSL	473443	100
 <p>Standard bore extension set with non-removable slide clamp, ULTRASITE valve injection site, SPIN-LOCK® connector. Components are not made with DEHP or natural rubber latex. Priming volume: 1.5 mL, Length: 8 in. (20 cm)</p>	US9060	474580	50
 <p>Standard bore extension set with two ULTRASITE valve injection sites 28 in. and 6 in. above distal end, non-removable slide clamp, SPIN-LOCK connector. Components are not made with DEHP or natural rubber latex. Priming volume: 5.5 mL, Length: 35 in. (89 cm)</p>	CSE33SL	473436	50