










IV ADMINISTRATION SETS

Description	Product Code	REF Number	Unit Per cs
INJECTION CAPS			
 <p>Intermittent Injection Cap - opaque yellow body, proximal self-sealing injection port, distal male luer lock connector with recessed luer taper. Components are not made with DEHP or natural rubber latex. Priming volume: 0.16 mL. Length: 0.875 in. (2.2 cm)</p>	IN1000	654237	1000
 <p>Intermittent Injection Cap - proximal injection port, clear body, distal male luer lock connector. Components are not made with DEHP or natural rubber latex. Priming volume: 0.2 mL. Length: 1.125 in. (2.9 cm)</p>	IN2000	418020	100
 <p>Intermittent Injection Cap - proximal injection port, distal luer slip connector. Components are not made with DEHP or natural rubber latex. Priming volume: 0.2 mL. Length: 1.125 in. (2.9 cm)</p>	IN3000	418025	100
 <p>Intermittent Injection Cap - proximal injection port, distal male luer lock connector. Components are not made with DEHP or natural rubber latex. Priming volume: 0.1 mL. Length: 0.75 in. (1.9 cm)</p>	IN4000	418030	100
REPLACEMENT CAPS			
 <p>Blue Replacement Cap - male luer lock.</p>	B1000	418017	100
Blue Replacement Cap - bulk pack of 1,000 in dispenser boxes of 100 each.	B1000B	418027	1000
White Replacement Cap - male luer lock.	C1000	418015	100
 <p>Replacement Caps - blue male luer lock cap and white female luer lock cap.</p>	BW1000	BW1000	100
 <p>White Replacement Cap - female luer lock.</p>	W1000	474900	100
 <p>RED CAP™ luer cap - Dual function red cap with male and female end. Bulk pack of 1,000 in dispenser</p>	R2000B	R2000B	1000
 <p>Blue Cap luer cap - Dual function blue cap with male and female end. Bulk pack of 1,000 in dispenser boxes of 100 each.</p>	B2000B	B2000B	1000