




## Blood Sets

### Y-Type Blood / Solution Sets with Standard Blood Filter (170 to 260 Micron) (10 drops / mL)

catalogue number	description	approx overall length	units / case
JC7751	Y-Type Blood / Solution Set with standard blood filter (170 to 260 micron), with male Luer lock adapter	98" (2.5 m)	48
			
Set has tubing that is compatible with the COLLEAGUE series, FLO-GARD 6201 and 6301 series volumetric infusion pumps			
JC6392	Y-Type Blood / Solution Set with standard blood filter (170 to 260 micron), 1 injection site, male Luer lock adapter	84" (2.1 m)	48
			
Set has tubing that is compatible with the COLLEAGUE series, FLO-GARD 6201 and 6301 series volumetric infusion pumps			
2C6750H	Y-Type Blood / Solution Set with standard blood filter (170 to 260 micron), 1 injection site, male Luer lock adapter	115" (2.9 m)	48
			
Set has tubing that is compatible with the COLLEAGUE series, FLO-GARD 6201 and 6301 series, SIGMA Spectrum volumetric infusion pumps			
JC6723	Y-Type Blood/Solution Set with large standard blood filter (170 to 260 micron), pressure pump, 2 injection sites, male Luer lock adapter	103" (2.6 m)	48
			
Set has tubing that is compatible with the COLLEAGUE series, FLO-GARD 6201 and 6301 series volumetric infusion pumps			