



Bellavita Auto Bath Lifter



FIGURE A



FIGURE B



HAND CONTROL



412100200

- The lightest bathlift on the market at only 20.5 lbs
- The Bellavita lowers to 2.3" (Figure A)
- Seat height reaches 18.8" for deeper tubs
- The backrest reclines to 50° at its lowest position
- Completely tool-free set up makes it easy-to-install and operate When not in use, back easily folds down for space-saving storage (Figure B)
- Watertight hand control floats
- The easy-to-use quick charge lithium-ion battery is accommodated in the hand control
- Two-piece construction makes the lift easy to assemble/disassemble for transport
- Padded seat and backrest provide added comfort
- Washable hygienic covers are sold separately
- Hygienic cutaway and slimline profile for easy handling
- 4 quick-release suction cups for anti-slip bath base
- Headrest available as an option
- Simplify transfers, positioning and rotating when using your bath lift with a Vitaturn Swivel and Transfer Aid
- Suction cups adhere to smooth surfaces only, and not textured surfaces

PRODUCTS

PRODUCTS	Description
477200252	Auto Bath Lifter w/white cover set, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Overall Dimensions	32.6" (L) (with backrest in a reclined position) x 27.5" (W) (seat incl. side flaps)
Seat Dimensions	13.7" (W) x 16.5" (D) x 22.4" (L)
Sitting Height	2.3" - 18.8"
Backrest Dimensions	25.1" (H) x 12.5" - 14.1" (W)
Backrest Reclines	50°
Baseplate Dimensions	12.2" (W) x 22" (L)
Load Capacity Max.	300 lbs.
Total Weight (excl. hand control)	20.5 lbs.
Shipping Weight	25.3 lbs.
Warranty	Limited 3 Year

ACCESSORIES

ACCESSORIES	Description
412100200	Rotating & Transfer Aid, White, 1/ea
410200252	Headrest - White, White, 1/ea
410200312	Headrest - Blue, Blue, 1/ea
460900205	Side Flap, 1/ea

PARTS

PARTS	Description
460900252	Cover-Set Comfort - White, White, 1/bx

PARTS	Description
460900312	Cover-Set Comfort - Blue, Blue, 1/bx
460900403	Charger for Bellavita, 1/bx
415252	Plug + Washer for Cover, White, 1/cv
415312	Plug + Washer for Cover, Blue, 1/cv
415432	Plug + Washer for Cover, Grey, 1/ea
460900213	Suction Cup with clip, Grey, 1/ea
460900220	Set of Suction Cups with Clips, Grey, 1/bx
460900600	Hand Control Including Storage Battery, 1/ea
460900225	Swiveling Sucker Set for Bellavita, Clear, 4/pkg
460900600	Hand Control Including Storage Battery, 1/ea
460900700	Back Assembly Complete, 1/ea
460900800	Seat Assembly Complete, 1/ea
460900403	Charger for Bellavita, 1/bx
460900209	Scissor Roller, 1/ea
460900213	Suction Cup with clip, 1/ea
460900220	Set of Suction Cups with Clips, 1/bx
460900225	Swiveling Sucker Set for Bellavita, 4/pkg
460900200	Scissors Complete, 1/ea
460900205	Side Flap, 1/ea
460900106	Motor, 1/ea
460900103	Gear Rack, 1/ea
460900223	Base Plate, 1/ea
460900211	Red Stop Lever, 1/ea
460900111	Two Wire Cable, 1/ea
410200252	Headrest - White, 1/ea
410200312	Headrest - Blue, 1/ea
410200432	Headrest - Gray, 1/ea