

Suction Catheter Trays & Solutions LATEX FREE



Suction Catheter Trays with Solution

- Sterile, coil packed suction catheter tray with chimney valve
- Two (2) walled blue nitrile latex-free exam gloves

CODE	FRENCH SIZE	TIP	SHIP CASE
------	-------------	-----	-----------

Sterile Saline – 100 mL Bottles

12062	6 Fr, graduated	DeLee	24
12082	8 Fr, graduated	DeLee	24
12102	10 Fr	Whistle	24
12142	14 Fr	Whistle	24
12153	14 Fr with drape	Whistle	24
12182	18 Fr	Whistle	24

Sterile Water – 100 mL Bottles

10062	6 Fr, graduated	DeLee	24
10082	8 Fr, graduated	DeLee	24
10102	10 Fr	Whistle	24
10122	12 Fr	Whistle	24
10142	14 Fr	Whistle	24
12191	14 Fr with drape	Whistle	24
10182	18 Fr	Whistle	24



Sterile Irrigating Solutions

- Sterile water and saline for suction and irrigation
- Safety seal

CODE	DESCRIPTION	SHIP CASE
------	-------------	-----------

Sterile Water Bottles - (8) Six Packs

10000	3.38 oz (100 mL) water	48
-------	------------------------	----

Sterile Saline Bottles - (8) Six Packs

1020	3.38 oz (100 mL) saline	48
------	-------------------------	----

Sterile Water Bottles - Bulk

1024	3.38 oz (100 mL) water	48
------	------------------------	----

Sterile Saline Bottles - Bulk

1022	3.38 oz (100 mL) saline	48
------	-------------------------	----



Mucus Traps



Argyle™ Mucus Traps with Vacuum Breaker and Filter

- For oral or wall/mechanical suction
- Advanced hydrophobic filter
- Whistle tip suction catheter with Slide-Text™ finish
- Transparent 20 cc graduated mucus trap with removable cap
- Additional cap for transport to lab
- Sterile

CODE	TRAP SIZE	CATHETER FRENCH SIZE	SHIP CASE
8888257485	20 cc	6 1/2 (2.2 mm)	50
8888257477	20 cc	8 (2.7 mm)	50
8888257469	20 cc	10 (3.3 mm)	50



Caution: This product contains natural rubber latex which may cause allergic reactions.



Argyle™ Mucus Traps with Vacuum Breaker

- For wall/mechanical suction
- Directional valve
- Whistle tip suction catheter with Slide-Text finish
- Transparent 20 cc graduated mucus trap with removable cap
- Additional cap for transport to lab
- Sterile
- Latex free

CODE	TRAP SIZE	CATHETER FRENCH SIZE	SHIP CASE
8888257386	20 cc	6 1/2 (2.2 mm)	50
8888257378	20 cc	8 (2.7 mm)	50
8888257360	20 cc	10 (3.3 mm)	50



Argyle™ Suction Catheter with Mucus Trap

- For oral suction
- Whistle tip catheter with Slide-Text finish
- Transparent 20 cc graduated mucus trap with removable cap
- Suction tube with molded mouthpiece
- Sterile
- Latex free

CODE	TRAP SIZE	CATHETER FRENCH SIZE	SHIP CASE
8888257550	20 cc	6 1/2 (2.2 mm)	50
8888257535	20 cc	8 (2.7 mm)	50
8888257527	20 cc	10 (3.3 mm)	50