

# CORFLO<sup>®</sup> NASOENTERIC

Opaque polyurethane feeding tube (unless otherwise noted), medication/flushing port, pre-inserted stylet, tungsten weighted tip, C-19<sup>™</sup> hydrophilic lubricant on tube tip and internal lumen, radiopaque

## CORFLO<sup>®</sup>-ULTRA PILL NG—WITH STYLET

Anti-clog, pill shape bolus, 3 gm\* adult weighted, with stylet

Code	NEW ENFit <sup>®</sup> Code	Description	Packaging
20-6438	40-6438	8 FR, 109cm Length	10 each
20-6431	40-6431	10 FR, 109cm Length	10 each

\*Approximate weight

## CORFLO<sup>®</sup>-ULTRA NG—WITH STYLET

Anti-clog, smooth shape bolus, 3 gm\* adult weighted, with stylet

Code	NEW ENFit <sup>®</sup> Code	Description	Packaging
20-7366	40-7366	6 FR, 91cm Length, 1.5 gm*	10 each
20-7368	40-7368	8 FR, 91cm Length	10 each
20-7361	40-7361	10 FR, 91cm Length	10 each
20-7438	40-7438	8 FR, 109cm Length	10 each
20-7431	40-7431	10 FR, 109cm Length	10 each
20-7432	40-7432	12 FR, 109cm Length	10 each

\*Approximate weight


## CORFLO<sup>®</sup>-ULTRA NG

Anti-clog, smooth shape bolus, 3 gm\* adult weighted, without stylet


Code	NEW ENFit <sup>®</sup> Code	Description	Packaging
20-5361	40-5361	10 FR, 91cm Length	10 each
20-5362	40-5362	12 FR, 91cm Length	10 each
20-5438	40-5438	8 FR, 109cm Length	10 each
20-5431	40-5431	10 FR, 109cm Length	10 each
20-5432	40-5432	12 FR, 109cm Length	10 each

\*Approximate weight


### ULTRA DESIGN



### SMOOTH DESIGN



### PILL DESIGN



NOTE: Please see package inserts for directions for use for each of these products.