

DeVilbiss iGo® Portable Oxygen Concentrator

HCPCS: E1390 + E1392



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Integrated PulseDose- oxygen conserving technology recognized for clinically effective pulse-dose oxygen delivery extends the use time of the DeVilbiss iGo Portable Oxygen System, offering a clinically sound and highly efficient ambulatory oxygen system Built-in OSD- (oxygen sensing device) ensures accurate oxygen delivery and reduced periodic maintenance schedule
- High-performance capabilities, compact design weighing less than 20 lbs with battery
- Fast battery charge time – 2 to 4.5 hours dependent upon the flow setting
- Easy-to-read battery life gauge
- OxyTrack Software provides an integrated solution for viewing performance and usage information on any DeVilbiss iGo Portable Oxygen System
- The external iGo Battery Charger is designed to charge batteries separate from the iGo unit (Figure A)
- Can be easily transported using the Deluxe Carry Bag (Figure C) or Wheeled Cart (Figure B)
- Optional accessory bag also available (Figure D-F)

PRODUCTS	Description
306DS	DeVilbiss iGo® Portable Oxygen Concentrator, EA
306DS-C	DeVilbiss iGo® Portable Oxygen Concentrator, EA

SPECIFICATIONS	Description
Altitude	0 -13,123 feet (0 - 4,000 meters) Tested @ ~933 hPa
Battery Warranty	1 Year Limited
Continuous Flow	1 to 3 LPM
Control Unit Airflow	8LPM
Control Unit Alarms	Visual/ Audible
Control Unit Cycle Time	10, 15, 20, 25 Minutes
Control Unit Type	Compressor
Mattress Cover Material	Fluid-Resistant Nylon
Mattress Dimensions	36" (W) X 80" (L) X 8" (D)
Max. Recommended Continuous Flow	3 LPM
Operating Temperature	41°F to 104°F (5°C to 40°C)
Oxygen Concentration	91% ± 3% For All Flow Settings
PulseDose® Settings	1 to 6
Sound Level	40 dBA
Unit Dimensions	15" (W) x 11" (L) x 8" (H)
Warranty	1 Year Limited
Weight Capacity (13046 13046-HD)	350 lbs.
Weight with battery	19 lbs. (8.6 kg)

SPECIFICATIONS	Description
Weight without battery	15.5 lbs. (7.0 kg)
iGo® Warranty	3 Year Limited

ACCESSORIES	Description
MC44D-605	EXT LIFE INTAKE BACT FILTR-DFT, EA
PV5LD-617	HOUR METER - DFT, EA
444-554	4 FT HOSE 1/8, EA
444-549	4 FT REINFORCED HOSE 1/2, EA
306DS-616	INTERNAL INTAKE FILTER, EA
PV5LD-651	FINAL BACTERIA FILTER - DFT, EA
306DS-652	Adapter, DC TO DC, EA

PARTS	Description
306CH	iGo Battery Charger, EA
306D-413	iGo Battery Pack, EA
306DS-601	115V CORD, EA
306DS-611	Gross particle filter 6/PK, EA
306DS-627	HUMIDIFIER KIT, EA
306DS-635	DELUXE IGO CARRY CASE, EA
306DS-651	AC TO DC ADAPTER, EA
306DS-655	iGo Accessory Bag , S. S. P., EA
4650D-609	PUSH-IN FEET (4-PK), PKG
R217P62	MAXO2+A OXYGEN ANALYZER, EA
306DS-602	230V EUROPE CORD, EA
306DS-603	230V UK CORD, EA
306DS-604	230V AU CORD, EA
306DS-605	230V CHINESE CORD, EA
306DS-606	CARTON WITH INSERTS, EA
444-507	ELBOW HUMID ADAPTER - DFT, EA