

### Flexam® Sterile Nitrile

#### Powder-Free Sterile Nitrile Examination Gloves

Flexam® Sterile Nitrile examination gloves offer all the comfort and protection of Stretchy Nitrile gloves, but are sterilized. They're soft and flexible, don easily and feel natural while in use. But they also provide the high level of protection you need in demanding clinical situations. They deliver freedom of movement and tactile sensitivity. Flexam® Sterile Nitrile exam gloves are available in either singles or pairs.

#### Features and Benefits:

- Stretchy nitrile formulation conforms to hands
- Microtextured surface for a sure grip
- 11 in. beaded cuff provides added protection
- Appropriate for those with natural rubber latex sensitivity
- Tested for use with various chemicals and chemotherapy drugs<sup>1</sup>

### Flexam® Sterile Nitrile



- Powder-free
- Not made with natural rubber latex
- Ambidextrous
- Color: Blue
- Nominal fingertip thickness: 5.9 mil (0.15mm)
- Nominal length: 11.8 in. (300mm)

Cat. No.	Size	Qty.
<b>Singles</b>		
N8820	Small	50/bx, 250/cs
N8821	Medium	50/bx, 250/cs
N8822	Large	50/bx, 250/cs
N8823	X-large	50/bx, 250/cs
<b>Pairs</b>		
N8830	Small	40/bx, 200/cs
N8831	Medium	40/bx, 200/cs
N8832	Large	40/bx, 200/cs
N8833	X-large	40/bx, 200/cs

<sup>1</sup> **Caution:** Gloves used for protection against chemotherapy drug and laboratory chemical exposure should be selected specifically for the type of chemicals being used. Note: Review Safety Data Sheets (SDS) for the chemicals/drugs being used to determine required level of protection.

### Triflex® Sterile Vinyl

#### Powder-Free Sterile Vinyl Examination Gloves

Triflex® Sterile Vinyl examination gloves are made from the same material as the Stretchy Synthetic gloves, offering tactile sensitivity and freedom of movement. Packaged in singles or pairs, the gloves are suitable for sterile tasks of short duration (less than 15 minutes) such as minor suture removals and dressing changes.

#### Features and Benefits:

- Vinyl formulation helps provide additional elasticity and fit
- Smooth finish
- Beaded cuff helps facilitate donning and helps reduce roll down
- Appropriate for those with natural rubber latex sensitivity

### Triflex® Sterile Vinyl

- Powder-free
- Not made with natural rubber latex
- Ambidextrous
- Color: Cream
- Nominal fingertip thickness: 4.7 mil (0.12mm)
- Nominal length: 9.5 in. (245mm)

Cat. No.	Size	Qty.
<b>Singles</b>		
2D7001PF	Small	100/bx, 200/cs
2D7002PF	Medium	100/bx, 200/cs
2D7003PF	Large	100/bx, 200/cs
<b>Pairs</b>		
2D7011PF	Small	50/bx, 200/cs
2D7012PF	Medium	50/bx, 200/cs
2D7013PF	Large	50/bx, 200/cs