

## Folding Steel Commode Retail



RETAIL BOX

- Easy, tool-free assembly
- Easily opens and folds for storage and transportation
- Durable plastic snap-on seat and lid
- Powder-coated steel-welded construction
- Plastic arms for added comfort
- Includes 7.5 qt commode bucket with carry handle, cover and splash shield
- Retail Packaged

### PRODUCTS

PRODUCTS	Description
RTL11158KDR	Folding Steel Commode, E0163, 1/ea

### Description

### SPECIFICATIONS

SPECIFICATIONS	Description
Width	18.125" between arms   22.25" Outside
Seat Dimensions	13.75" (W) x 16.25" (D)
Seat Height	17.0" – 23.0"
Weight	13 lbs.
Carton Shipping Weight	17.5 lbs.
Weight Limit	350 lbs.
Warranty	Limited Lifetime

### Description