


Titanium Skin Protection and Positioning Gel/Foam Wheelchair Cushion



FPT-1



FPT-1

- The Titanium cushion is specifically designed to offer superior pressure redistribution and positioning
- High-resilient, high-density foam combined with a dual chamber gel bladder encourages proper patient posture through immersion
- Middle contour acts as a built-in leg abductor
- Low-shear, fluid-resistant stretch cover helps to maintain sensitive skin tissue
- Non-skid bottom offers additional stability
- Meets flammability standard TB117

PRODUCTS

Description

PRODUCTS	Description
FPT-1	3" Gel/Foam Cushion, E2607, 1/ea
FPT-2	3" Gel/Foam Cushion, E2607, 1/ea
FPT-1818	3" Gel/Foam Cushion, E2607, 1/ea
FPT-2018	3" Gel/Foam, E2607, 1/ea
FPT-3	3" Gel/Foam Cushion, E2607, 1/ea
FPT-4	3" Gel/Foam Cushion, E2608, 1/ea
FPT-5	3" Gel/Foam Cushion, E2608, 1/ea
FPT-6	3" Gel/Foam Cushion, E2608, 1/ea

SPECIFICATIONS

Description

SPECIFICATIONS	Description
Cover Material	Fluid Resistant Stretch Top Non-Skid Bottom
Weight Capacity	Up to 20" 275 lbs. 22" and Up 500 lbs.
Warranty	1 Year Limited