

COMPRIFOAM®

FOAM BANDAGE



Comprifoam® open cell foam bandage is designed to evenly distribute the pressure from compression bandages over extremities, thus preventing constrictions and ensuring an effective compression therapy.

INDICATIONS

- Padding under compression bandages

FEATURES AND BENEFITS

- Breathable
- Non-sliding
- Tear resistant
- Easy for patients to self bandage
- Anatomically conforming allowing freedom of movement
- Resistant to body oils, lotions and perspiration
- Reusable
- May be steam sterilized

Dimensions	Code	Box Contents	Selling Unit
10 cm x 2.5 m x 0.4 cm	7529400	1 roll	Box
12 cm x 2.5 m x 0.4 cm	7529500	1 roll	Box

ARTIFLEX®

NON-WOVEN PADDING BANDAGE



Artiflex® is a multipurpose padding bandage that ensures even distribution of pressure in phlebological or lymphedema applications. Also useful as under padding of plaster or synthetic casts.

INDICATIONS

- Padding under compression bandages for sensitive skin
- Under padding of plaster and synthetic casts

FEATURES AND BENEFITS

- Non-absorbent
- Applies easily without creasing
- Provides added comfort in wrapping for lymphedema
- Breathable
- Autoclavable

Dimensions	Code	Box Contents	Selling Unit
10 cm x 3 m	09046	30 rolls	Box
15 cm x 3 m	09047	20 rolls	Box

FOAM RUBBER PAD



Dimensions	Code	Box Contents	Selling Unit
9 cm x 1 cm	78496	Unit	Box