

Nasogastric Feeding Tubes — Adult



Kangaroo™ Feeding Tubes with Dobbhoff™ Tip

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert[†]
- Purple color to denote enteral feeding
- Weighted feeding tubes with and without flow-through stylets
- Radiopaque polyurethane material
- CM markings
- MRI conditional once stylet is removed
- DEHP-Free

| Code | Includes Stylet | Rigid Outlet Port | Weight | French Size | Length | SHIP Case |
|------------|-----------------|-------------------|--------|-------------|--------------|-----------|
| 8884710867 | Yes | Yes | 5g | 8 | 43" (109 cm) | 10 |
| 8884710842 | No | Yes | 7g | 8 | 43" (109 cm) | 10 |
| 8884710859 | Yes | Yes | 7g | 8 | 43" (109 cm) | 10 |
| 8884710826 | Yes | Yes | 7g | 8 | 55" (140 cm) | 10 |
| 8884711006 | Yes | Yes | 7g | 10 | 43" (109 cm) | 10 |
| 8884711246 | No | No | 7g | 12 | 43" (109 cm) | 10 |
| 8884711253 | Yes | Yes | 7g | 12 | 43" (109 cm) | 10 |