



Battery-Powered Patient Lift with Rechargeable, Removable Battery



FIGURE A



FIGURE B




FIGURE C

- Same features as 13240
- Includes removable, rechargeable battery, item #13240BATT (Figure A) external charger item #13255 (Figure B) and wall-mounted-bracket item #13260 (Figure C)
- Ideal for institutional use
- To avoid downtime – order an additional rechargeable battery, Item #13240BATT

PRODUCTS

PRODUCTS	Description
13242	Battery-Powered Patient Lift, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Boom Height	27.5" - 76"
Base Width	41" (open) 24" (closed)
Base Height	4.5"
Base Clearance	2.25"
Casters	4.5" Rear
Weight	105 lbs.
Carton Shipping Weight	120 lbs.
Weight Limit	450 lbs.

PARTS

PARTS	Description
13255	External Charger, 1/ea
13240BATT	Rechargeable Battery, 1/ea
13260	Wall-Mounted Bracket, 1/ea
13050	4/6 point cradle, 1/cs
13240-CC	Cradle Clip (6 ea), 6/pkg
13240CAB	Cradle Attachment Bolt, 1/ea
13240FP	Frame Plug, 1/ea
13240MU	Upper Mast Boom, 1/ea
13240UBA	Upper Bolt For Actuator, 1/ea
13240MHB	Mast Hinge Bolt, 1/ea
13240BATT	Rechargeable Battery, 1/ea
13240BC	Battery Cell * *Need to order 2 each, 1/ea
13240PH	Push Handle, 1/ea
13255	External Charger, 1/ea
13240A	Actuator, 1/ea
13240HC	Hand Control, 1/ea
13240BSB	Spreader Bar, 1/ea
13240SBH	Spreader Bar Handle, 1/ea

PARTS	Description
13240M	Mast, 1/ea
13240F	Footplate, 1/ea
13240FGP	Foot Pad, 1/ea
13240CAL	Caster with Lock, 1/ea
13240BAR-PR	Base Actuating Rod, 1/pr
13240B	Base, 1/cv
13240CA	Caster without Lock, 1/ea
13240PC	Power Cord for 13240 charger, 1/ea
13255PC	Power Cord for 13255 external charger, 1/ea
13240RPA	Receiver Plate Assm., 1/ea
13240ARB	Receiver Plate Bolt, 1/pkg
13240ARB-N	Ser # 4A1206 & Later, 1/pkg
13058PB	Cradle Pivot Bolt, 1/ea
13260	Wall-Mounted Bracket, 1/ea
13240-CC-SL	Cradle Clips Spring-Loaded
1324614	Control Box
1324626	Charger Stand
13255PCH	Power Cord Hanger
CPCFUSE	Fuse