

# Patient Room Sharps Collectors

BD Patient Room Collectors are available in BD Recykleen™ brand and BD standard brand and come in a variety of sizes and colors. All BD Patient Room Collectors have visible fill-level indicators and handles for safe transportation.

## Key product features include:

- Compatible with locking cabinets to provide security and prevent unwanted access to collector contents
- Counterbalanced door facilitates hands-free disposal and is designed to reduce the potential of sharps injuries
- Designed to improve tamper resistance by preventing hand access to the contents
- Available in BD Recykleen™ brand



305053



305517



## Patient / Exam Room

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
<b>BD Recykleen™ Sharps Collectors</b>				
305053	3 gal BD Recykleen™ Counterbalanced Door Pearl	16.6 x 10.7 x 6	10	305097, 305099, 305409, 305588, 305552
305054	3 gal BD Recykleen™ Counterbalanced Door Red	16.6 x 10.7 x 6	10	
305517	5.4 qt Counterbalanced Door Red	12 x 12 x 4.8	20	305096, 305098, 305447, 305021
305551	5.4 qt Counterbalanced Door Clear	12 x 12 x 4.8	20	
305426	5.4 qt Counterbalanced Door Red	12 x 12 x 4.8	12	
305427	5.4 qt Counterbalanced Door Clear	12 x 12 x 4.8	12	
305422	2 gal Counterbalanced Door Pearl	12.45 x 10.72 x 6	10	305097, 305099, 305409, 305588, 305552
305435	2 gal Counterbalanced Door Red	12.45 x 10.72 x 6	10	
305436	3 gal Counterbalanced Door Red	16.6 x 10.7 x 6	10	
305443	5.4 qt Counterbalanced Side Entry Red	11.7 x 11.6 x 4.5	20	305447, 305021
305444	5.4 qt Counterbalanced Side Entry Pearl	11.7 x 11.6 x 4.5	20	
305425	5.4 qt Counterbalanced Side Entry Pearl	11.7 x 11.6 x 4.5	12	