


PrimeMat 2.0 Impact Reduction Fall Mats



PM20BA

- Full 72" length provides maximum coverage alongside bed
- Tough and lightweight - Mats are constructed from high quality molded Polyurethane foam inner core (210kgs/cubic meter) to provide maximum absorption and rebound to reduce impact in the event of a fall
- Passed 5 foot "Egg Drop Test" - an egg was dropped from 5 feet above the mat without damage
- Low profile with tapered edges allow equipment such as overbed tables or patient lifts to roll over easily
- Anti-slip bottom netting keeps the mat secure to help prevent accidental tripping
- Embossed non-slip PVC leather look top is not only attractive but extremely durable and puncture resistant
- Non-slip pattern helps minimize slipping and is easy to clean with mild soap and water
- Available in two colors; Brown & Gray

PRODUCTS

Description

PM20BA	Fall Mat 2.0 - Brown, 1/ea, 5/cs
PM20GA	Fall Mat 2.0 - Gray, 1/ea, 5/cs

SPECIFICATIONS

Description

Unit Dimensions	24" (W) x 72" (L) x 0.75" (H)
Warranty	1 Year Limited