

ProBasics™ Rollators



A. ProBasics™ Flame Finish Aluminum Rollator with Loop Brakes

Featuring a lightweight, aluminum welded frame, full swivel front casters and easy squeeze brakes for progressive braking control, this rollator is easy for individuals to maneuver. The flip-up, padded seat reveals a quick-access storage pouch and the entire rollator can be easily folded for storage and transportation. HCPCS Code: E0143, E0156

Model #	Description
1037	Aluminum Rollator - Flame Red
1037BL	Aluminum Rollator - Flame Blue
1037BK	Aluminum Rollator - Flame Silver
1037F	Aluminum Rollator - Flame Fuchsia



B. ProBasics™ Deluxe Aluminum Rollator

Ideal for those who spend more time outdoors, the Deluxe Aluminum Rollator features a durable, lightweight aluminum frame, large 8" wheels and a rugged in-line braking system. This rollator offers a seamless, padded polyurethane seat for improved comfort and ease of care. The textured seat is waterproof, anti-bacterial and flips up to reveal a large, zippered compartment for secure storage. HCPCS Code: E0143, E0156

Model #	Description
1032MAR	Deluxe Aluminum Rollator - Flame Red
1032BLMAR	Deluxe Aluminum Rollator - Flame Blue



C. ProBasics™ Bariatric Aluminum Rollator

The bariatric aluminum rollator features a durable welded frame, locking loop brakes, large 8" wheels and a 400 lb. weight capacity to provide individuals with a heavy duty, reliable product. With a low seat-to-floor height and 5 ribs on the heavy duty front and rear wheels, individuals can feel more stable. The Bariatric Aluminum Rollator folds easily for transportation and storage and comes with a large storage pouch. HCPCS Code: E0149, E0156

Model #	Description
1033MAR	Bariatric Aluminum Rollator - Flame Red
1033MARBL	Bariatric Aluminum Rollator - Flame Blue

Model #	Dimensions inside Hand Grips		Seat Dimensions	Seat-to-Floor Height	Handle Height Adjustment	Caster Dimensions	Product Weight	Weight Capacity	Limited Warranty
1037, 1037BL 1037BK 1037F	18"	14"W 12"D	21"	31"-37"	6"	14 lb.	300 lb.	3 Years Frame 1 Year Seat	
1032MAR 1032BLMAR	18"	14"W 12"D	21"	31.5"-35.5"	8"	16 lb.	300 lb.	3 Years Frame 1 Year Seat	
1033MAR 1033MARBL	23.5"	19"W 13"D	22"	33.5"-38.5"	8"	19.4 lb.	400 lb.	3 Years Frame 1 Year Seat	

NOTE: Rollators are not intended to be used as wheelchairs

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.