

Suction Cup Grab Bars



RETAIL BOX



INSET B



INSET C



INSET D



- Installs and removes without tools or professional installation
- Can be installed at any angle
- Large suction cups provide an extremely strong hold; suction cups are 3.75" in diameter
- Release levers make installing and removing the suction cup grab bars quick and easy
- Red/Green color indicators on top of the handle shows "safe secure" suction or to "re-affix" suction
- Avoid grout lines for proper placement
- Will not damage property
- Use to assist with balance only
- Do not use to leverage body weight
- Always check to make sure the bar is securely attached before using
- Uninstall and reinstall as needed
- Retail Packaged

PRODUCTS

Description

RTL13083	White & Beige 12" Suction Cup Grab Bar, 3/cs
RTL13085	Black & Chrome 12" Suction Cup Grab Bar, 3/cs

SPECIFICATIONS

Description

Length	12"
Weight	.5 lbs.
Carton Shipping Weight	3.5 lbs.
Warranty	Limited Lifetime