

AirLife cushion nasal cannulas



AirLife cushion cannula 002600

These nasal cannulas have a soft, clear cushion design that supports the nasal prongs to optimize patient comfort. Cushion cannulas are available in sizes from neonate to adult with a variety of oxygen tubing lengths. All oxygen supply tubing is crush-resistant to minimize occlusion and kinking.

Cat. no.	Description	Qty.
Neonate, infant and pediatric cannulas		
002611	Neonate, with 7' (2.1m) tubing	50/cs
002601	Infant, with 7' (2.1m) tubing	50/cs
002602	Pediatric, with 7' (2.1m) tubing	50/cs
Adult cannulas		
002600	With 7' (2.1m) tubing	50/cs

Vital Signs Capnoflex end tidal CO₂ sampling cannulas



End tidal CO₂ sampling cannula for Capnoflex

The Vital Signs ETCO₂ sampling cannulas for Capnoflex have a proprietary design for use exclusively with the Capnoflex module utilized with the GE Dash 3000/4000/5000™ patient monitors. They are offered in a variety of styles for all patients, in both nasal and oral/nasal designs with or without oxygen delivery. Many of the Capnoflex kits have Nafion sampling lines offering extended usage for patients that are receiving high humidity. All cannulas come with 8' (2.4m) sampling lines.

Cat. no.	Description	Qty.
2013066-003	Infant CO ₂ sampling nasal cannula only	10/cs
2067767-003	Infant CO ₂ sampling nasal cannula with dehumidification sampling line	10/cs
2013066-007	Infant CO ₂ sampling nasal cannula with O ₂ delivery	10/cs
2067767-006	Infant CO ₂ sampling nasal cannula with dehumidification sampling line and O ₂ delivery	10/cs
2013066-002	Pediatric CO ₂ sampling nasal cannula only	10/cs
2067767-002	Pediatric CO ₂ sampling nasal cannula with dehumidification sampling line	10/cs
2013066-005	Pediatric CO ₂ sampling nasal cannula with O ₂ delivery	10/cs
2067767-005	Pediatric CO ₂ sampling nasal cannula with dehumidification sampling line and O ₂ delivery	10/cs
2013067-002	Pediatric CO ₂ sampling oral/nasal cannula only	10/cs
2067767-008	Pediatric CO ₂ sampling oral/nasal cannula with dehumidification sampling line	10/cs
2013067-004	Pediatric CO ₂ sampling oral/nasal cannula with O ₂ delivery	10/cs
2067767-010	Pediatric CO ₂ sampling oral/nasal cannula with dehumidification sampling line and O ₂ delivery	10/cs
2013066-001	Adult CO ₂ sampling nasal cannula only	10/cs
2067767-001	Adult CO ₂ sampling nasal cannula with dehumidification sampling line	10/cs
2013066-004	Adult CO ₂ sampling nasal cannula with O ₂ delivery	10/cs
2067767-004	Adult CO ₂ sampling nasal cannula with dehumidification sampling line and O ₂ delivery	10/cs
2013067-001	Adult CO ₂ sampling oral/nasal cannula only	10/cs
2067767-007	Adult CO ₂ sampling oral/nasal cannula with dehumidification sampling line	10/cs
2013067-003	Adult CO ₂ sampling oral/nasal cannula with O ₂ delivery	10/cs
2067767-009	Adult CO ₂ sampling oral/nasal cannula with dehumidification sampling line and O ₂ delivery	10/cs
2013069-001	CO ₂ sampling line with male luer connection - 6.5 ft (2m)	10/cs
2013069-002	CO ₂ dehumidification sampling line with male luer connection - 6.5 ft (2m)	10/cs

What are your quality goals?

- Improving patient satisfaction
- Improving patient safety
- Standardizing products and process
- Reducing cost

If you checked any, or all, of the above contact your sales representative today to learn how AirLife® can help you achieve your quality objectives.