

## LUMBROSACRAL SUPPORTS

**PROCARE®****ComfortFORM™ Back**

Overlapping double pull straps help support the lumbar and abdominal area to promote proper spinal alignment and help prevent and protect injuries including strains and sprains. Silicone gripper strips on either side of back pad help prevent migration. Optional shoulder straps available. Latex free.

**Indication:** Helps prevent and protect injuries including strains and sprains; Aids in proper alignment while sitting and standing for prolonged periods. Support may be worn under or over clothing.

	Waist Measurement	Size
79-89352	20" – 25"	XS
79-89353	25" – 30"	S
79-89355	30" – 34"	M
79-89357	34" – 38"	L
79-89358	38" – 42"	XL
79-89359	42" – 56"	2XL
79-89359-10	56" – 62"	3XL
79-89350	Shoulder Straps	Univ.





**PROCARE**

**Sacro-Lumbar Support w/Compression Straps**

Low profile 9" elastic support with double side pull compression straps. Neoprene posterior panel accommodates ComfortFORM™ Moldable Insert #79-89120. Foam lumbar pad is included with 79-8250x style. No stays.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

Black	White	Waist Meas.	Size
79-82502		18" – 24"	XS
79-82503	79-89263	24" – 33"	S
79-82505	79-89265	33" – 39"	M
79-82507	79-89267	39" – 45"	L
79-82508	79-89268	45" – 53"	XL
79-82509	79-89269	53" – 59"	2XL
79-82509-10		59" – 65"	3XL
79-82509-11		65" – 72"	4XL



**PROCARE**

**Clinic Retention Support with Compression Straps**

Foam/pile, with 9" panel elastic construction and double side pull elastic compression straps. Two sizes. Accommodates ComfortFORM™ moldable insert #79-89120.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89240	24" – 40"	S/M
79-89241	40" – 54"	L/XL



**PROCARE**

**10" Criss-Cross Support**

Universal style criss-cross elastic support with two plastic posterior stays.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89040	28" – 50"	Univ.



**PROCARE**

**Criss-Cross Support with Compression Straps**

10" elastic criss-cross support with compression straps. Four 9" stays—two aluminum that can be easily formed to help ensure proper support and two plastic to maintain position of support.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89183	24" – 30"	S
79-89185	30" – 36"	M
79-89187	36" – 42"	L
79-89188	42" – 48"	XL
79-89189	48" – 52"	2XL



**PROCARE**

**10" Double-Pull Sacro-Lumbar Support**

Premium all-elastic construction with compression straps. Four 9" posterior plastic stays for stability and support.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89003	31" – 34"	S
79-89005	35" – 38"	M
79-89006	36" – 42"	M/L
79-89007	39" – 42"	L
79-89008	43" – 46"	XL
79-89009	47" – 54"	2XL
79-89009-10	55" – 60"	3XL



**PROCARE**

**Sacro-Lumbar Support**

Elastic construction with felt back pad. Four posterior preformed aluminum stays. Three adjustable web compression straps with front clip buckle closure for custom fit.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

13"	Waist Measurement	Size
79-89023	28" – 32"	S
79-89024	32" – 36"	M
79-89026	36" – 40"	M/L
79-89027	40" – 44"	L
79-89028	44" – 48"	XL
79-89029	48" – 52"	2XL



**PROCARE**

**Sacro-Lumbar Support with Mesh Back**

10" elastic construction with ventilated mesh back panel. Four 9" plastic posterior stays. Side pull elastic compression straps with contact closure.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89203	24" – 30"	S
79-89205	30" – 36"	M
79-89207	36" – 42"	L
79-89208	42" – 48"	XL
79-89209	48" – 54"	2XL



**PROCARE**

**Neoprene/Elastic Sacro-Lumbar Support w/Compression Straps**

Low profile 10" elastic construction with neoprene back panel and side pulls for warmth, compression and adjustability. Two plastic side stays prevent rolling and 3 spiral posterior stays provide added support. Available in sized and universal.

**Indication:** Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89013	28" – 32"	S
79-89015	34" – 38"	M
79-89016	36" – 40"	M/L
79-89017	40" – 44"	L
79-89018	44" – 48"	XL
79-89019	48" – 52"	2XL
79-89010	32" – 52"	Univ.