




Round-Handle Canes with Foam Grip



- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring
- Locking ring prevents rattling (Figure A)
- Handle height adjusts from 28.75" to 37.75"
- 10345C-6 Contains 2 each Silver (Figure B), Bronze (Figure C) and Black (Figure D)
- Weight Capacity 300 lbs

PRODUCTS	Description	
10345C-6	Contains 2 Each Silver, Black & Bronze., 6/cs	
RTL10342	Silver, 6/cs	
SPECIFICATIONS	Description	
Handle Height	28.75" - 37.75"	
Weight Limit	300 lbs.	
Warranty	Limited Lifetime	