



2020 Product Catalog

ALL ENFit® ALL the TIME!





Vesco

Latin verb meaning "to feed"

Feeding Innovation

Our promise to deliver a continuous stream of break-through
enternal products to better serve our customers





High-Quality Tube-Feeding Devices

- ISO 9001 & 13485 Certified Quality Management System
- FDA Registered & 510(k) Premarket Clearance
- Made without natural rubber latex, DEHP or other phthalates plasticizers



Unbeatable Service

- Friendly, skilled and experienced customer care specialists
- Full portfolio of ENFit products direct or through distributors
- Responsive when and where you need us



Competitive Pricing

- Lean, agile organization driven to keep costs down
- High-quality products without the premium price
- Volume discounts for even more savings




ENFit Tip Syringes

ENFit Tip Syringes for Feed, Flush & Irrigation (5, 10, 20, 35 & 60mL)


- ISO 80369-3 compliant ENFit connector designed to reduce the risk of tubing misconnections
- Enteral specific syringes for use in initial tube placement, flushing/hydration, bolus feeding and irrigation
- Printed graduated length markings in mL
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers
- One cap included in each blister pack

ENFit Low Dose Tip Syringes (1 & 3 mL)


- Designed to minimize process changes in pharmacy and bedside administration
- Demonstrated dose accuracy administering smaller volumes compliant with GEDSA*
- Internal lumen consistent with traditional male syringe design reduces deadspace and priming volume
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers
- Sterile blister packs with one cap included in each pack



ENFit Standard Tip



ENFit Low Dose Tip (LDT)



ENFit Syringe with Transition Connector

THUMB RING IRRIGATION ENFIT SYRINGES	Item #	Units/Carton
60mL	VED-661	180
60mL + 500mL Container	VED-665	180

ENFit is a federally registered trademark of GEDSA in multiple jurisdictions throughout the world. Used with permission.

* GEDSA August, 2016 Guidance Document

Syringe Item Numbers


ENFIT TIP SYRINGES			CLEAN (Non-Sterile)		
Size (mL)	STERILE		with Cap	with Transition Connector	Units/Carton
	Item #	Units/Case			
1 (LDT)	VED-601EO w/ Caps	800	N/A	N/A	N/A
3 (LDT)	VED-603EO w/ Caps	800	N/A	N/A	N/A
5	VED-605EO	800	VED-605	VED-605TC	500
10	VED-610EO	800	VED-610	VED-610TC	500
20	N/A	N/A	VED-620	VED-620TC	500
35	VED-635EO	400	VED-635	VED-635TC	500
60	VED-660EO	400	VED-660	VED-660TC	180



ENFit Pharmacy Items

ENFit Medication Bottle Fill Caps (Multiple Sizes)

- Pharmacy friendly purple cap designed to aid in single handed placement
- Various sizes to meet most common medication bottles
- Fits inside bottle opening to allow continued use of child proof caps
- Ideal for aseptic technique
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers



ENFit Transition Connectors

- Transition Connector consistent with the design and taper of a catheter tip syringe
- Available prepackaged within the inner case or within each syringe package
- Also available as a stand alone item by the case

Transition Connectors

ENFit Syringe Tip Cap (Fits All Sizes)

- Pharmacy friendly purple cap designed to aid in single handed placement
- Ideal for aseptic technique and secure for transport
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers



ENFit Medicine Straw (50 mm)

- Connects to an ENFit Tip Syringe to easily facilitate the withdrawal of medication or feed

Important Safety Information

- Use ENFit tip syringes only with ENFit feeding tube ports. Feeding tubes with different open female ports are not compatible with Vesco ENFit tip syringes
- ENFit tip syringes have the potential to misconnect with small bore connectors of other healthcare applications
- Single-patient use and Rx only



ENFit Medicine Straw

Pharmacy Item Numbers


ENFit Medication Bottle Adapters			
Size #	Inner Diameter (mm)	Item #	Units/Case
1	12 - 15	VED-690	100
2	14.5 - 16.5	VED-691	100
3	17 - 20.5	VED-692	100
4	19 - 23	VED-693	100
5	24 - 26	VED-694	100
OTHER			Units/ Carton
Transition Connectors			VED-001 300
ENFit Medicine Straw			VED-670 100
ENFit Syringe Screw-on Cap			VED-697 4000
ENFit Syringe Pop-on Cap			VED-699P 4000
ENFit Syringe Pop-on Cap			VED-699PF 4000



VED-697



VED-699P



VED-699PF

ENFit SYRINGE PHARMACY BULK PACK				
Size (mL)	Item #	Units/ Case	AMBER Item #	Units/ Case
1 (LDT)	VED-601B	4000	VED-601AB	4000
3 (LDT)	VED-603B	4000	VED-603AB	4000
5	VED-605B	2000	VED-605AB	1000
10	VED-610B	1000	VED-610AB	1000
20	N/A	N/A	VED-620AB	1000
35	VED-635B	800	VED-635AB	800
60	VED-660B	400	VED-660AB	400

Caps sold separately for sterile/non-sterile bulk



Low Profile Extension Sets

Universal Low Profile Feeding Tube Extension Sets


- Bolus Feed Single ENFit Port
- Continuous Feed Dual ENFit Port

Features

- Compatible with MIC Key® and MINI One® low profile gastrostomy feeding tubes
- Dual ENFit Ports with 90 degree connector allows for continuous feeding, easy flushing and administration of liquid medication without disconnecting set
- Single ENFit port & straight connector for bolus feeding, flushing and medication delivery
- Not made with natural rubber latex, Bisphenol A or other phthalate plasticizers

Important Safety Information

- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes
- Device is non sterile
- See individual product labeling for directions and important safety information
- Extension Set has the potential to misconnect with small bore connectors of other healthcare applications
- Single patient use and Rx only



Extension Set Item Numbers

UNIVERSAL EXTENSION SETS Set Description	12 Inches (30cm)	24 Inches (60cm)	Units/Case
	Item #	Item #	
CONTINUOUS Dual Port Right Angle	VED-1000	VED-1001	10
BOLUS Single Port Straight	VED-1002	VED-1003	10
BOLUS Single Port Right Angle	VED-1004	VED-1005	10



Gravity Feeding Sets

Gravity Feeding Screw Cap and Bag Sets with ENFit connector and Transition Connector


Features

- Gravity Bag Sets
 - 3 Times Stronger and up to 35% Faster than competing PVC bags*
 - Snap cap with rigid neck for easy filling and closure
 - Easy-to-read graduations and easy view translucent material
- Closed System Sets
 - Screw cap allows for flexibility in feeding containers with a 40mm opening
 - Convenient suspension bag allows for use with smaller containers (baby bottles) for easy hanging
- Not made with natural rubber latex, DEHP, Bisphenol A or other phthalate plasticizers
- Exit port allows complete product drainage
- Device intended to deliver liquid nutrition formulas or water to a patient's feeding tube

ENFit Transition Set

- Pre-attached ENFit Transition Connector allows compatibility with existing feeding tube female ports
- Removing ENFit Transition Connector allows compatibility with new ISO standard connection on new feeding tubes with ENFit connectors

* Performance testing vs. competing PVC bags based on engineering studies by Vesco Medical Prod. Dev.



1000mL Gravity Feed Set

Gravity Set Item Numbers

GRAVITY FEED SETS	Item #	Units/Case
Screw Cap ENFit w/ Transition Connectors	VED-020	24
ENPlus Safety Screw w/ Transition Connectors	VED-026	30
1000mL ENFit w/ Transition Connectors	VED-032	30
500mL ENFit w/ Transition Connectors	VED-038	30
1000mL TOP FILL (PVC) w/ Transition Connectors	VED-045	30
1000mL TOP FILL (PVC)	VED-046	30

Important Safety Information

- To avoid risk of a tubing misconnection, **ONLY** connect the distal connector to an ENTERAL feeding tube. This device is compatible with a flexible female connector on the patient's enteral feeding tube. Not for IV use
- This device should be replaced at least every 24 hours
- Do not store for prolonged periods at extreme temperatures
- Device is non-sterile
- See individual product labeling for directions and important safety information
- Single patient use and Rx only




Gastrostomy Feeding Tubes

Dual ENFit Port Replacement Gastrostomy Feeding Tubes (14 to 24 Fr)

Features

- Allows for easy flushing and administration of liquid medication without disconnecting set
- Medical-grade silicone tubing with clear shaft and graduation marks in 1cm increments
- Rounded skin disk for retention
- Gastric balloon for internal retention with rounded tip, open distal end and radiopaque tip marker
- Not made with natural rubber latex, Bisphenol A or other phthalate plasticizers



Dual ENFit
Port G-Tube

Single ENFit
Port G-Tube

Single ENFit Port Replacement Gastrostomy Feeding Tubes (16 to 22 Fr)

Features

- Biocompatible-grade silicone tubing with clear shaft and graduation marks in 2cm increments.
- Rounded skin disk with bumpers for retention
- Gastric balloon for internal retention with rounded flush tip and open distal end
- Not made with natural rubber latex, Bisphenol A or other phthalate plasticizers

Important Safety Information

- Use only administration/feeding sets with ENFit connectors
- Device is sterilized (EO for Ethyl-o-xylene)
- Tube has the potential to misconnect with small bore connectors of other healthcare applications
- Single patient use and Rx only

G-Tube Item Numbers

G TUBES	Single Port			Dual Port			
	French Size	Item #	Balloon Size (mL)	Units/Case	Item #	Balloon Size (mL)	Units/Case
14	N/A	N/A	N/A	N/A	VED-214	5	5
16	VED-116	10	2	VED-216	5	5	
18	VED-118	15	2	VED-218	20	5	
20	VED-120	15	2	VED-220	20	5	
22	VED-122	15	2	VED-222	20	5	
24	N/A	N/A	N/A	VED-224	20	5	



Nasoenteric Feeding Tubes

Single and Dual ENFit Port Nasoenteric Feeding Tube with and without Stylet


- French Sizes: 5, 6, 8, 10, 12, 14 & 16
- Length Options: 18 (45cm), 36 inches (91 cm) & 45 inches (114cm)

Features

- Single and Dual Port ENFit connectors for feeding and administering medication
- Premium medical-grade kink resistant polyurethane tubing
- 4 side ports and distal end port improves flow and reduces likelihood of tube clogging
- Water-activated lubricant aids in tube placement and allows for aseptic lubrication
- Radiopaque for X-ray confirmation of tube placement
- Not made with natural rubber latex, DEHP, Bisphenol A or other phthalate plasticizers
- Centimeter graduated length markings
- Available with and without flow through stylet for ease of placement and confirmation

Important Safety Information


- Use only administration/feeding sets with ENFit connectors. Feeding sets with different connectors are not compatible with Vesco feeding tubes
- Device is non sterile
- Tube has the potential to disconnect with small bore connectors of other healthcare applications
- Single patient use and Rx only



Single & Dual ENFit Port NG Tubes

NE Tube Item Numbers

NE TUBES	36 Inches (91 cm)				45 Inches (114 cm)				Units/Case
	Single Port		Dual Port		Dual Port		Dual Port		
	No Stylet	Stylet	No Stylet	Stylet	No Stylet	Stylet	No Stylet	Stylet	
French Size	Item #		Item #		Item #		Item #		Item #
5	VED-309	N/A	N/A	N/A	N/A	N/A	N/A	N/A	10
6	VED-304	N/A	VED-306	VED-306S	VED-406	VED-406S	VED-406	VED-406S	10
8	VED-300	VED-300S	VED-308	VED-308S	VED-408	VED-408S	VED-408	VED-408S	10
10	VED-301	VED-301S	VED-310	VED-310S	VED-410	VED-410S	VED-410	VED-410S	10
12	VED-302	VED-302S	VED-312	VED-312S	VED-412	VED-412S	VED-412	VED-412S	10
14	VED-303	VED-303S	VED-314	VED-314S	VED-414	VED-414S	VED-414	VED-414S	10
16	N/A	N/A	VED-316	VED-316S	N/A	N/A	N/A	N/A	10



Vesco Medical

Tube feeding experts delivering high quality products with unbeatable service at a competitive price. Our portfolio is All ENFit All The Time!

Group Purchasing Agreements:

HealthTrust (HT31810)
Premier (PP-DI-1468)

MHA (MS60238N)
Vizient (MS6216)

Distributor Partners:

Concordance
Cardinal/Indemed

McKesson
IMS

Owens Minor
Medline

Preferred Medical
Rightway

www.VescoMedical.com
info@vescomedical.com
(855) 239-0131

GEDSA

ENFit is a federally registered trademark of GEDSA in multiple jurisdictions throughout the world. Used with permission.
©2020 Vesco Medical, LLC MKT 1102-8-08



For more information regarding new ISO Standard ENFit connectors, visit www.stayconnected.org



1039 Kingsmill Parkway
Columbus, OH 43229