



Open Suction
1-800-962-9888



OPEN SUCTION

Table of Contents

Open Suction Products

Tips, valves and packaging1

Single Suction Catheters

Single Chimney Catheters2

Single Directional Catheters 2

Sterile Suction Catheter Trays Without Solutions

Suction Catheter Kits with Chimney Valve3

Chimney Suction Catheter Trays3

Directional Suction Catheter Trays.....3

Sterile Suction Catheter Kits, Trays & Solutions

Chimney Suction Catheter Kits
with Sterile Saline (Bottle)4

Chimney Suction Catheter Trays
with Sterile Water (Bottle)4

Sterile Water Bottles4

Sterile Saline Bottles.....4

Mucus Traps

Argyle™ Mucus Traps5

Specimen Containers

Argyle™ Lukens Specimen Container6

Argyle™ Specimen Trap.....6

Specimen Collection Set.....6

LUKI™ Tubes6

Tracheostomy Care

Tracheostomy Care Trays7

Tracheostomy Tube Holder.....7

Index

Index 8

Open Suction Products

Covidien suction catheters are available in two tip styles: the smooth, atraumatic DeLee tip with staggered eyes, and the angled Whistle tip. Chimney valves and directional valves deliver easy fingertip control with virtually no residual vacuum in the open position.

TIPS

DeLee

- Standard tip for all Covidien pediatric (5 Fr – 8 Fr) suction catheters
- Smooth tip minimizes trauma

Whistle

- Round, open tip provides efficient non-traumatic flow rate


VALVES

Chimney Valve

- Minimal residual vacuum in the open position

Directional Valve

- Minimizes risk of cross contamination



SUCTION CATHETER PACKAGING

Straight



Coiled



Single Suction Catheters LATEX FREE



Argyle™ Single Suction Catheters with Chimney Valve

- Chimney valves utilize Whistle tip for adult sizes (10-18 Fr) and DeLee tip for Pediatric (5-8 Fr)
- Sterile, single-use, disposable
- 1 cm graduation marks start at 4 cm on pediatric and 16 cm on adult catheters

CODE	FRENCH SIZE	TIP	SHIP CASE
------	-------------	-----	-----------

Straight Packed

30500	5 Fr pediatric, graduated	DeLee	50
30688	6 Fr pediatric, graduated	DeLee	50
30888	8 Fr pediatric, graduated	DeLee	50
31000	10 Fr	Whistle	50
31044	10 Fr, graduated	Whistle	50
31200	12 Fr	Whistle	50
31244	12 Fr, graduated	Whistle	50
31400	14 Fr	Whistle	50
31444	14 Fr, graduated	Whistle	50
31445	14 Fr, with Coudé tip	Whistle	50
31600	16 Fr	Whistle	50
31644	16 Fr graduated	Whistle	50
31800	18 Fr	Whistle	50

Coil Packed

30520	5 Fr pediatric, graduated	DeLee	50
30620	6 Fr pediatric, graduated	DeLee	50
30820	8 Fr pediatric, graduated	DeLee	50
31020	10 Fr	Whistle	50
31220	12 Fr	Whistle	50
31420	14 Fr	Whistle	50
31620	16 Fr	Whistle	50
31820	18 Fr	Whistle	50

Straight Packed with Straight Connector (without Control Valve)

33425	14 Fr	Whistle	50
-------	-------	---------	----

Coil Packed with Straight Connector (without Control Valve)

33000	10 Fr	Whistle	50
33400	14 Fr	Whistle	50
33800	18 Fr	Whistle	50

Argyle™ Single Suction Catheters with Directional Valve

- Directional valve has ISO standard colors which easily identify French sizes
- Sterile, single-use, disposable
- 1 cm graduation marks start at 4 cm on pediatric and 16 cm on adult catheters

CODE	FRENCH SIZE	COLOR	TIP	SHIP CASE
------	-------------	-------	-----	-----------

Straight Packed

141906	6 Fr pediatric, graduated	Light Green	DeLee	100
141908	8 Fr pediatric, graduated	Light Blue	DeLee	100
141900	10 Fr	Black	DeLee	100
141901	12 Fr	White	DeLee	100
141902	14 Fr	Green	DeLee	100
141903	16 Fr	Orange	DeLee	100
141904	18 Fr	Red	DeLee	100



Suction Catheter Trays LATEX FREE



31279



37024



150982

Argyle™ Suction Catheter Mini-Soft-Kits without Solution

- Sterile, coil packed suction catheter with chimney valve
- Two (2) walled blue nitrile latex-free exam gloves
- Pop-up solution cup in sterile inner wrap

CODE	FRENCH SIZE	TIP	SHIP CASE
------	-------------	-----	-----------

30579	5 Fr, graduated	DeLee	50
30677	6 Fr, graduated	DeLee	50
30877	8 Fr, graduated	DeLee	50
31077	10 Fr, graduated	DeLee	50
31079	10 Fr	Whistle	50
31279	12 Fr	Whistle	50
31479	14 Fr	Whistle	50
31679	16 Fr	Whistle	50
31879	18 Fr	Whistle	50

Argyle™ Suction Catheter Tray with Chimney Valve

- Sterile suction catheter kit with chimney valve
- Two (2) Walled blue nitrile latex-free exam gloves

CODE	FRENCH SIZE	TIP	SHIP CASE
------	-------------	-----	-----------

Coil Packed

36524	5 Fr, graduated	DeLee	50
36626	6 Fr, graduated	DeLee	50
36826	8 Fr, graduated	DeLee	50
37024	10 Fr	Whistle	50
37224	12 Fr	Whistle	50
37424	14 Fr	Whistle	50
37524	16 Fr	Whistle	50
37724	18 Fr	Whistle	50

Argyle™ Suction Catheter Tray with Directional Valve

- Sterile, coil packed suction catheter with directional valve
- Two (2) walled blue nitrile latex-free exam gloves
- Mini tray offers maximum efficiency

CODE	FRENCH SIZE	TIP	SHIP CASE
------	-------------	-----	-----------

140988	8 Fr, pediatric, graduated	DeLee	100
140980	10 Fr	DeLee	100
140981	12 Fr	DeLee	100
140982	14 Fr	DeLee	100
140983	16 Fr	DeLee	100

Suction Catheter Trays & Solutions LATEX FREE



Suction Catheter Trays with Solution

- Sterile, coil packed suction catheter tray with chimney valve
- Two (2) walled blue nitrile latex-free exam gloves

CODE	FRENCH SIZE	TIP	SHIP CASE
------	-------------	-----	-----------

Sterile Saline – 100 mL Bottles

12062	6 Fr, graduated	DeLee	24
12082	8 Fr, graduated	DeLee	24
12102	10 Fr	Whistle	24
12142	14 Fr	Whistle	24
12153	14 Fr with drape	Whistle	24
12182	18 Fr	Whistle	24

Sterile Water – 100 mL Bottles

10062	6 Fr, graduated	DeLee	24
10082	8 Fr, graduated	DeLee	24
10102	10 Fr	Whistle	24
10122	12 Fr	Whistle	24
10142	14 Fr	Whistle	24
12191	14 Fr with drape	Whistle	24
10182	18 Fr	Whistle	24



Sterile Irrigating Solutions

- Sterile water and saline for suction and irrigation
- Safety seal

CODE	DESCRIPTION	SHIP CASE
------	-------------	-----------

Sterile Water Bottles - (8) Six Packs

10000	3.38 oz (100 mL) water	48
-------	------------------------	----

Sterile Saline Bottles - (8) Six Packs

1020	3.38 oz (100 mL) saline	48
------	-------------------------	----

Sterile Water Bottles - Bulk

1024	3.38 oz (100 mL) water	48
------	------------------------	----

Sterile Saline Bottles - Bulk

1022	3.38 oz (100 mL) saline	48
------	-------------------------	----



Mucus Traps



Argyle™ Mucus Traps with Vacuum Breaker and Filter

- For oral or wall/mechanical suction
- Advanced hydrophobic filter
- Whistle tip suction catheter with Slide-Text™ finish
- Transparent 20 cc graduated mucus trap with removable cap
- Additional cap for transport to lab
- Sterile

CODE	TRAP SIZE	CATHETER FRENCH SIZE	SHIP CASE
8888257485	20 cc	6 1/2 (2.2 mm)	50
8888257477	20 cc	8 (2.7 mm)	50
8888257469	20 cc	10 (3.3 mm)	50



Argyle™ Mucus Traps with Vacuum Breaker

- For wall/mechanical suction
- Directional valve
- Whistle tip suction catheter with Slide-Text finish
- Transparent 20 cc graduated mucus trap with removable cap
- Additional cap for transport to lab
- Sterile
- Latex free

CODE	TRAP SIZE	CATHETER FRENCH SIZE	SHIP CASE
8888257386	20 cc	6 1/2 (2.2 mm)	50
8888257378	20 cc	8 (2.7 mm)	50
8888257360	20 cc	10 (3.3 mm)	50



Argyle™ Suction Catheter with Mucus Trap

- For oral suction
- Whistle tip catheter with Slide-Text finish
- Transparent 20 cc graduated mucus trap with removable cap
- Suction tube with molded mouthpiece
- Sterile
- Latex free

CODE	TRAP SIZE	CATHETER FRENCH SIZE	SHIP CASE
8888257550	20 cc	6 1/2 (2.2 mm)	50
8888257535	20 cc	8 (2.7 mm)	50
8888257527	20 cc	10 (3.3 mm)	50



Specimen Containers



8888258608

Argyle™ Lukens Specimen Container

- 14 Fr Whistle tip suction catheter with Slide-Tex finish
- 20 cc specimen trap graduated at 2 cc intervals
- Extra trap cap
- Suction tube with integral tapered connector
- Patient label
- Sterile
- Latex free

CODE	CATHETER FRENCH SIZE	SHIP CASE
8888258608	14 (4.7 mm)	50



8884724500

Argyle™ Specimen Traps

- Screw-on cap
- Removable male connector
- Full 40 cc or 70 cc capacity with 1 cc graduations
- Patient label
- Sterile
- Latex free

CODE	CAPACITY	SHIP CASE
8884724500	40 cc	50
8884724504	70 cc	50

LATEX Caution: This product contains natural rubber latex which may cause allergic reactions.



15000

Specimen Collection Set

- Save-T-Vac™ valve
- Screw on cap
- Patient label
- Sterile

CODE	CAPACITY	SHIP CASE
15000	30 mL	25

Luki™ Disposable Aspirating Tube

LATEX Caution: This product contains natural rubber latex which may cause allergic reactions.

- Sterile
- Disposable

CODE	CAPACITY	TUBE SIZE	SHIP CASE
8886864304	8 cc	3.25"	48
8886864604	20 cc	6.25"	48



Tracheostomy Care LATEX FREE

Covidien provides hospitals and alternate care facilities with a wide selection of tracheostomy care trays which offer the finest in nursing convenience and infection control. Designed with the guidance of experienced nursing personnel, these superior total-care trays provide exactly what is required for efficient and quick rendering of trach care procedures.

Tracheostomy Care Trays

- Sterile, Single-use

Item Code	42201	47800	47802	47815	47835	47850	47885	47890	47891	47892
Ship Case	20	20	24	20	20	20	20	20	48	50
Tray	•	•	•	•	•	•	•		•	
Blue Nitrile Gloves (pair)	•	•	•	•	•	2	•			
Pipe Cleaners	3	2	2	2	2	2	2	2	2	2
Trach Brush	•	•	•	•	•	•	•	•		
Cotton Tip Applicators	3	2	2	2	2	2	2	2	2	2
Drape	•	•	•	•	•	•	•			
Twill Tape Trach Holder	•	•	•	•	•	•	•	•		
Removable Basin	3	•	•	•	•	•	•			
Fenestrated Trach Sponge	3	•	•	•	•	•	•		2	2
Suction Catheter (14 Fr)				•						
Plastic Forceps					•	•				
Gauze	4	4	4	4	4	4	4	2	2	2
Hydrogen Peroxide			2				2		•	•
Mask						•				
Sterile Saline			•			•				
CSR Wrap	•									



47802



47891



47892



Tracheostomy Tube Holders

CODE	DESCRIPTION	SHIP CASE
49000	Non-sterile	100

INDEX

CODE	PRODUCT	HCPCS CODES	PAGE	CODE	PRODUCT	HCPCS CODES	PAGE
10000	Sterile Water Bottles - (8) Six Packs	A9999	4	31444	Sterile Straight Packed Suction Catheters, Chimney Valve 14 Fr	A4624	2
10062	Suction Catheter Trays with Sterile Water 6 Fr	A4624, A4216, A4930	4	31445	Sterile Straight Packed Suction Catheters, Chimney Valve 14 Fr	A4624	2
10082	Suction Catheter Trays with Sterile Water 8 Fr	A4624, A4216, A4930	4	31479	Suction Catheter Mini Soft Kits without Solution 14 Fr	A4624, A4930	3
10102	Suction Catheter Trays with Sterile Water 10 Fr	A4624, A4216, A4930	4	31600	Sterile Straight Packed Suction Catheters, Chimney Valve 16 Fr	A4624	2
10122	Suction Catheter Trays with Sterile Water 12 Fr	A4624, A4216, A4930	4	31620	Sterile Single Coil Packed Suction Catheters, Chimney Valve 16 Fr	A4624	2
10142	Suction Catheter Trays with Sterile Water 14 Fr	A4624, A4216, A4930	4	31644	Sterile Straight Packed Suction Catheters, Chimney Valve 16 Fr	A4624	2
10182	Suction Catheter Trays with Sterile Water 18 Fr	A4624, A4216, A4930	4	31679	Suction Catheter Mini Soft Kits without Solution 16 Fr	A4624, A4930	3
1020	Sterile Saline Bottles - (8) Six Packs	A9999	4	31800	Sterile Straight Packed Suction Catheters, Chimney Valve 18 Fr	A4624	2
1022	Sterile Saline Bottles - Bulk	A9999	4	31820	Sterile Single Coil Packed Suction Catheters, Chimney Valve 18 Fr	A4624	2
1024	Sterile Water Bottles - Bulk	A9999	4	31879	Suction Catheter Mini Soft Kits without Solution 18 Fr	A4624, A4930	3
12062	Suction Catheter Trays with Sterile Saline 6 Fr	A4624, A4930, A4216	4	33000	Sterile Coil Packed Suction Catheters, Straight Connector 10 Fr	A4624	2
12082	Suction Catheter Traays with Sterile Saline 8 Fr	A4624, A4930, A4216	4	33400	Sterile Coil Packed Suction Catheters, Straight Connector 14 Fr	A4624	2
12102	Suction Catheter Trays with Sterile Saline 10 Fr	A4624, A4930, A4216	4	33425	Sterile Straight Packed Suction Catheters, Straight Connector 14 Fr	A4624	2
12142	Suction Catheter Trays with Sterile Saline 14 Fr	A4624, A4930, A4216	4	33800	Sterile Coil Packed Suction Catheters, Straight Connector 18 Fr	A4624	2
12153	Suction Catheter Trays with Sterile Saline 14 Fr	A4624, A4930, A4216	4	36524	Argyle Suction Catheter Tray with Chimney Valve 5 Fr	A4624	3
12182	Suction Catheter Trays with Sterile Saline 18 Fr	A4624, A4930, A4216	4	36626	Argyle Suction Catheter Tray with Chimney Valve 6 Fr	A4624	3
12191	Suction Catheter Trays, Sterile Water 14 Fr with drape	A4624, A4216, A4930	4	36826	Argyle Suction Catheter Tray with Chimney Valve 8 Fr	A4624	3
140980	Directional Suction Catheter Mini-Tray 10 Fr	A4624	3	37024	Argyle Suction Catheter Tray with Chimney Valve 10 Fr	A4624	3
140981	Directional Suction Catheter Mini-Tray 12 Fr	A4624	3	37224	Argyle Suction Catheter Tray with Chimney Valve12 Fr	A4624	3
140982	Directional Suction Catheter Mini-Tray 14 Fr	A4624	3	37424	Argyle Suction Catheter Tray with Chimney Valve 14 Fr	A4624	3
140983	Directional Suction Catheter Mini-Tray 16 Fr	A4624	3	37524	Argyle Suction Catheter Tray with Chimney Valve 16 Fr	A4624	3
140988	Directional Suction on Catheter Mini-Tray 8 Fr	A4624	3	37724	Argyle Suction Catheter Tray with Chimney Valve 18 Fr	A4624	3
141900	Sterile Straight Packed Suction Catheters, Directional Valve 10 Fr	A4624	2	42201	Tracheostomy Care Tray	A4629	7
141901	Sterile Straight Packed Suction Catheters, Directional Valve 12 Fr	A4624	2	47800	Tracheostomy Care Tray	A4629	7
141902	Sterile Straight Packed Suction Catheters, Directional Valve 14 Fr	A4624	2	47802	Tracheostomy Care Tray	A4629	7
141903	Sterile Straight Packed Suction Catheters, Directional Valve 16 Fr	A4624	2	47815	Tracheostomy Care Tray	A4629, A4624	7
141904	Sterile Straight Packed Suction Catheters, Directional Valve 18 Fr	A4624	2	47835	Tracheostomy Care Tray	A4629	7
141906	Sterile Straight Packed Suction Catheters, Directional Valve 6 Fr	A4624	2	47850	Tracheostomy Care Tray		7
141908	Sterile Straight Packed Suction Catheters, Directional Valve 8 Fr	A4624	2	47885	Tracheostomy Care Tray	A4629	7
15000	Specimen Collection		6	47890	Economy Trach Care Kit	A9270	7
30500	Sterile Straight Packed Suction Catheters, Chimney Valve 5 Fr	A4624	2	47891	Tracheostomy Care Tray	A9270	7
30520	Sterile Single Coil Packed Suction Catheters, Chimney Valve 5 Fr	A4624	2	47892	Soft Pack Trach Care Kit	A9270	7
30579	Suction Catheter Mini Soft Kits without Solution 5 Fr	A4624, A4930	3	49000	Tracheostomy Tube Holder	A4999	7
30620	Sterile Single Coil Packed Suction Catheters, Chimney Valve 6 Fr	A4624	2	8884724500	Argyle Specimen Traps	A9999	6
30677	Suction Catheter Mini Soft Kits without Solution 6 Fr	A4624, A4930	3	8884724504	Argyle Specimen Traps	A9999	6
30688	Sterile Straight Packed Suction Catheters, Chimney Valve 6 Fr	A4624	2	8886864304	LUKI Disposable Aspirating Tube		6
30820	Sterile Single Coil Packed Suction Catheters, Chimney Valve 8 Fr	A4624	2	8886864604	LUKI Disposable Aspirating Tube		6
30877	Suction Catheter Mini Soft Kits without Solution 8 Fr	A4624, A4930	3	8888257360	Argyle DeLee Mucus Traps with Vacuum Breaker	A9270	5
30888	Sterile Straight Packed Suction Catheters, Chimney Valve 8 Fr	A4624	2	8888257378	Argyle DeLee Mucus Traps with Vacuum Breaker	A9270	5
31000	Sterile Straight Packed Suction Catheters, Chimney Valve 10 Fr	A4624	2	8888257386	Argyle DeLee Mucus Traps with Vacuum Breaker	A9270	5
31020	Sterile Single Coil Packed Suction Catheters, Chimney Valve 10 Fr	A4624	2	8888257469	Argyle DeLee Mucus Traps with Vacuum Breaker and Filter	A9270	5
31044	Sterile Straight Packed Suction Catheters, Chimney Valve 10 Fr	A4624	2	8888257477	Argyle DeLee Mucus Traps with Vacuum Breaker and Filter	A9270	5
31077	Suction Catheter Mini Soft Kits without Solution 10 Fr	A4624, A4930	3	8888257485	Argyle DeLee Mucus Traps with Vacuum Breaker and Filter	A9270	5
31079	Suction Catheter Mini Soft Kits without Solution 10 Fr	A4624, A4930	3	8888257527	Argyle DeLee Suction Catheters with Mucus Trap	A9270	5
31200	Sterile Straight Packed Suction Catheters, Chimney Valve 12 Fr	A4624	2	8888257535	Argyle DeLee Suction Catheters with Mucus Trap	A9270	5
31220	Sterile Single Coil Packed Suction Catheters, Chimney Valve 12 Fr	A4624	2	8888257550	Argyle DeLee Suction Catheters with Mucus Trap	A9270	5
31244	Sterile Straight Packed Suction Catheters, Chimney Valve	A4624	2	8888258608	Argyle Lukens Specimen Container	A9270	6
31279	Suction Catheter Mini Soft Kits without Solution 12 Fr	A4624, A4930	3				
31400	Sterile Straight Packed Suction Catheters, Chimney Valve 14 Fr	A4624	2				
31420	Sterile Single Coil Packed Suction Catheters, Chimney Valve 14 Fr	A4624	2				

DISCLAIMER: The Centers for Medicare and Medicaid Services was the primary source of any HCPCS information contained in this document. Covidien does not represent or guarantee that any HCPCS information contained in this document is complete, accurate, or applicable to any particular patient or third-party payer and disclaims all liability for any consequence resulting from any reliance on such information. The final billing-related decisions for any product must be made by the health care provider considering the medical necessity of the product furnished as well as the requirements of third-party payers and any local, state, or federal laws and regulations that apply to the product furnished. Covidien is providing this information in an educational capacity with the understanding that Covidien is not engaged in rendering legal, accounting, or other professional services. Covidien encourages all health care providers to consult with their own advisors regarding coding and payment.

The products shown in this catalog may not be licensed or registered in all parts of the world.
Please contact your sales representative for availability

Cover photo credit: Thinkstock.

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.
Other brands are trademarks of a Covidien company. ©2012 Covidien.

H6200SUCT – 2M – 0612

