

Cutimed®

NEW DIMENSIONS IN INTEGRATED WOUND CARE



BSN medical delivers an integrated management pathway
with positive patient outcomes in all phases of complex wound healing.

Our innovative technologies offer a complete portfolio of patient and user-friendly products to effectively
remove bacteria, help rebuild tissue, and reduce pressure to successfully manage difficult wounds.

THERAPIES. HAND IN HAND.

www.bsnmedical.com

BSN medical

NEW DIMENSIONS IN INTEGRATED WOUND CARE

REMOVE Bacteria

with *Cutimed® Siltec Sorbact®*
featuring *DACC Technology*



- Quickly Binds and Inactivates Bacteria¹
- Improves Healing Time²
- Safe and Effective Even Against MRSA³

REBUILD Tissue

with help from *Cutimed® Epiona*
Collagen Dressing



- Promotes Tissue Granulation and Epithelialization⁴
- Native 90% Collagen/10% Alginate
- 97% Identical to Human Collagen⁴

REDUCE Pressure

with *Cutimed® Off-Loader*
Select TCC System



- The Gold Standard for TRUE Total Contact Casting
- Locks Heel and Toes in Place
- Achieves Greater Patient Compliance

Ordering Information

Ref-No.	Product	Size	Dressing/Box	Suggested HCPCS* Codes
7322700	CUTIMED® EPIONA	2" x 2" (5.1 x 5.1 cm)	10	A6021
7322701	CUTIMED® EPIONA	4" x 4" (10.2 x 10.2 cm)	10	A6021
7322702	CUTIMED® EPIONA	8" x 8" (20.3 x 20.3 cm)	10	A6023
7992900	CUTIMED® SILTEC SORBACT®**	3" x 3"	10	A6212
7992901	CUTIMED® SILTEC SORBACT®**	5" x 5"	10	A6212
7992902	CUTIMED® SILTEC SORBACT®**	6" x 6"	10	A6212
7992903	CUTIMED® SILTEC SORBACT®**	7" x 7"	5	A6213
7993300	CUTIMED® SORBACT® HYDROACTIVE B**	2.8" x 3.3"	10	A6242
7993302	CUTIMED® SORBACT® HYDROACTIVE B**	5.5" x 5.5"	10	A6242
7993303	CUTIMED® SORBACT® HYDROACTIVE B**	7.5" x 7.5"	10	A6242
7800901	CUTIMED® OFF-LOADER SELECT	COMPLETE KIT	1	CPT 29445
ASK REP	CANVAS ROCKER BOTTOM CAST SHOE	S, M, L, XL	1	TBD

1. Ljungh, A., Yanagisawa, N., and Wadstrom, T. Using the principle of hydrophobic interaction to bind and remove wound bacteria. *The Journal of Wound Care*, April 2006, Vol. 15, No. 4.
 2. Corsi, Alessandro. Comparison of a novel non-mediated bacteria-binding dressing to silver dressings in the management of acute and chronic skin lesions.
 3. Ronner, A.C., et al. Adhesion of methicillin-resistant Staphylococcus aureus to DACC-coated dressings. *Journal of Wound Care*, Vol. 23 / No. 10, October 2014.
 4. Data on file.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

* This document includes codes that might be used to bill for BSN medical products and related services. It is the provider's sole responsibility to determine and submit appropriate codes, charges, and modifiers for services rendered. BSN medical Inc. cannot guarantee coverage or reimbursement.

** Please consult your BSN sales representative for a complete list of Cutimed® Sorbact® sizes, dressings/box, and suggested HCPCS.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
 Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
 To order toll-free: BSN medical (+1) 800 552 1157

62325 R1 ©2015 BSN medical Inc. REV 10/15

