

Deluxe Aluminum Bath Chair With Tool-free Removable Back



INSET A



INSET B



RETAIL BOX

- Easy, tool-free assembly of back seat and legs (Figure A - B)
- Aluminum frame is lightweight, durable and corrosion proof
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability

PRODUCTS

PRODUCTS	Description
RTL12202KDR	Gray, 1/cs
12202KDRT-1	Teal, 1/cs
12202KDRY-1	Yellow, 1/cs
12202KDRB-1	Blue, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Outside Legs	19.5" (W) x 19.5" (D)
Seat Dimensions	19.25" (W) x 11.5" (D)
Seat Height	15.5" - 19.5"
Weight	7 lbs.
Carton Shipping Weight	8.5 lbs.
Weight Limit	400 lbs.
Warranty	Limited Lifetime