



DermaGinate™ Calcium Alginate Dressing

DermaGinate is a calcium alginate wound dressing for moderate to high exuding wounds. This highly absorbent dressing is derived from natural seaweed. Conforms readily to any wound bed. Forms a soothing, gel-like consistency on contact with moisture, yet maintains integrity for convenient removal. Assists in maintaining a moist wound healing environment.

***New!** In addition to our standard Rope dressing, we are pleased to offer our new DermaGinate 12" Fluffy Rope, which leaves less residue than standard calcium alginate dressings, and allows for easier filling of irregularly shaped wounds.*

Benefits:

- Absorbs exudate while minimizing the risk of maceration
- Supports moist wound healing & autolytic debridement
- Biocompatible
- Hemostatic
- Conforms easily
- Latex free

Indications:

A primary wound dressing for moderate to high exuding wounds, such as:

- pressure ulcers
- leg ulcers
- diabetic ulcers
- surgical wounds
- donor sites
- lacerations and abrasions

Product	Size	Reorder #	Units
DermaGinate	2x2	00260E	10/Box
DermaGinate	4x4	00270E	10/Box
DermaGinate	4x8	00272E	5/Box

Product (cont.)	Size	Reorder #	Units
DermaGinate	12" Rope	00275E	5/Box
DermaGinate	12" Fluffy Rope	00274E	5/Box

