

Multi-use Nestable Sharps Collectors

BD Multi-use Nestable Collectors are available in three sizes and can be used with accessories for added security and stability.

Key product features include:

- Clear tops with multiple opening styles to accommodate different sizes of waste and environments
- Handles designed for safe transport
- Funnel-style lids limit accessibility to the collector contents and minimize spillage if the collector tips over
- Open-top lids feature larger diameter openings to accommodate disposal needs in non-patient / public areas



305457

8 qt



Regular Funnel Clear Top Hinge Cap with Petals

8 qt



Open Clear Top Hinge Cap

14 qt and 6 gal



Large Funnel Clear Top Hinge Cap with Petals

14 qt



Open Clear Top with Tethered Cap



305343 with 305570 Stabilizer



305646 with 305485 Bracket



Multi-use Nestable

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305343	8 qt Hinge Cap	10.3 x 11.3 x 6.8	24	305485, 305018, 305570
305344	8 qt Hinge Cap with Petals	10.3 x 11.3 x 6.8	24	
305646	14 qt Tethered Cap	11.5 x 12.8 x 8.8	20	
305456	14 qt Hinge Cap	11.5 x 12.8 x 8.8	20	
305464	14 qt Hinge Cap with Needle Ports	11.5 x 12.8 x 8.8	20	
305480	14 qt Hinge Cap with Petals	11.5 x 12.8 x 8.8	20	
305457	6 gal Hinge Cap	17.5 x 12.8 x 8.8	12	
305465	6 gal Hinge Cap with Petals	17.5 x 12.8 x 8.8	12	
305078	6 gal Tethered Cap	17.5 x 12.8 x 8.8	12	

Brackets and Stabilizers

305485	Bracket	6 x 3.5 x 1.25	12	8 qt, 14 qt, 6 gal
305018	Wire Bracket	9.5 x 10.25 x 7.75	6	All 8 qt
305570	Stabilizer	2 x 10 x 7	5	All 8 qt
305494	Security Lock	1.75 x 1 x .25	12	305485, 305018, 305570

See Accessories on page 14.