

Pinniped Pediatric Commode

HCPCS: E1399



INSET A



PC 1000 R

- All-in-one welded steel frame available in red and blue
- Plastic seat and contoured back
- Durable snap-on toilet seat
- Height adjustable seat allows child to firmly place their feet on the floor for comfort, safety and security
- Adjustable hook & loop chest strap for additional support
- Molded plastic armrests provide additional comfort
- Rubber non-skid tips ensure that the legs will not mark the floor
- Complete with 7.5 qt removable commode bucket with metal handle and cover

PRODUCTS

Description

PC 1000 R	Red, 1/bx
PC 1000 BL	Blue, 1/bx

SPECIFICATIONS

Description

Seat Depth	15"
Seat Width	13"
Seat Height	10.5" – 14.5"
Armrest Height (from floor)	15.5" - 19.5"
Armrest Pad	2" (W) x 9" (D)
Overall Width	17"
Toilet Opening	6.5" (W) x 8" (D)
Weight	9.5 lbs
Carton Shipping Weight	13 lbs
Weight Limit	250 lbs
Warranty	Limited Lifetime

ACCESSORIES

Description

RTL12085	Biodegradable Sanitary Commode Liner, 12/bx, 6 bx/cs
----------	--